



MODEL:  
**POC-W22A-H81**

**21.5" Medical Panel PC with 4<sup>th</sup> Generation Intel® Dual-Core CPU, Touchscreen, Anti-bacteria Cover, Dual USB 3.0, Dual GbE LAN, RS-232/422/485, HDMI, VGA, Front IP 65 Rating and RoHS**

# User Manual

# Revision

---

Date	Version	Changes
January 24, 2017	1.01	Minor update
October 17, 2016	1.00	Initial release

# Copyright

---

## COPYRIGHT NOTICE

The information in this document is subject to change without prior notice in order to improve reliability, design and function and does not represent a commitment on the part of the manufacturer.

In no event will the manufacturer be liable for direct, indirect, special, incidental, or consequential damages arising out of the use or inability to use the product or documentation, even if advised of the possibility of such damages.

This document contains proprietary information protected by copyright. All rights are reserved. No part of this manual may be reproduced by any mechanical, electronic, or other means in any form without prior written permission of the manufacturer.

## TRADEMARKS

All registered trademarks and product names mentioned herein are used for identification purposes only and may be trademarks and/or registered trademarks of their respective owners.

## CONTACT INFORMATION

### **IEI Integration Corp.**

**Address:** No. 29, Zongxing Rd., Xizhi Dist.,  
New Taipei City 221, Taiwan

**Phone:** +886-2-8691-6798

**Fax:** +886-2-6616-0028

**Web Site:** [www.ieiworld.com](http://www.ieiworld.com)

**Sales Email:** [sales@ieiworld.com.tw](mailto:sales@ieiworld.com.tw)

# Manual Conventions

---



## **WARNING**

Warnings appear where overlooked details may cause damage to the equipment or result in personal injury. Warnings should be taken seriously.



## **CAUTION**

Cautionary messages should be heeded to help reduce the chance of losing data or damaging the product.



## **NOTE**

These messages inform the reader of essential but non-critical information. These messages should be read carefully as any directions or instructions contained therein can help avoid making mistakes.



## **HOT SURFACE**

This symbol indicates a hot surface that should not be touched without taking care.



## **OPERATING INSTRUCTION**

Follow operating instructions or consult instructions for use.



## **IEC 60417-5009: STAND-BY**

# Table of Contents

<b>1 INTRODUCTION.....</b>	<b>1</b>
1.1 OVERVIEW.....	2
1.2 MODEL VARIATIONS .....	3
1.3 FEATURES.....	4
1.4 FRONT PANEL.....	5
1.4.1 Backlit Touch Buttons .....	6
1.5 SIDE PANELS .....	7
1.6 BOTTOM PANEL.....	7
1.7 REAR PANEL.....	8
1.8 SYSTEM SPECIFICATIONS.....	9
1.9 DIMENSIONS.....	13
<b>2 UNPACKING .....</b>	<b>14</b>
2.1 UNPACKING.....	15
2.2 PACKING LIST.....	16
2.3 OPTIONAL ITEMS.....	17
<b>3 INSTALLATION .....</b>	<b>20</b>
3.1 ANTI-STATIC PRECAUTIONS .....	21
3.2 INSTALLATION PRECAUTIONS .....	21
3.3 INSTALLATION AND CONFIGURATION STEPS .....	22
3.4 REMOVING THE BACK COVER .....	22
3.5 HDD INSTALLATION.....	24
3.6 UPS BATTERY .....	26
3.7 HANDSET INSTALLATION (OPTIONAL) .....	27
3.7.1 Using VoIP Handset.....	28
3.8 HANDLE INSTALLATION (OPTIONAL).....	31
3.8.1 Barcode Reader Installation .....	32
3.8.2 Reading Light.....	34
3.9 3-IN-1 COMBO READER INSTALLATION (OPTIONAL) .....	35
3.10 USING RFID READER (OPTIONAL).....	37



3.11 RS-232/422/485 SERIAL PORT CONNECTION .....	40
3.12 AT/ATX MODE SELECTION .....	40
3.12.1 AT Power Mode.....	40
3.12.2 ATX Power Mode .....	41
3.13 CABLE COVER INSTALLATION .....	41
3.14 MOUNTING THE SYSTEM .....	42
3.14.1 Wall Mounting.....	43
3.14.2 Arm Mounting .....	46
3.14.3 Stand Mounting .....	47
3.15 POWERING ON THE SYSTEM.....	48
3.16 RESET THE SYSTEM .....	49
<b>4 DRIVER INSTALLATION.....</b>	<b>50</b>
4.1 AVAILABLE SOFTWARE DRIVERS .....	51
4.2 INTEL® CHIPSET DRIVER .....	52
4.3 INTEL® GRAPHICS DRIVER .....	52
4.4 AUDIO DRIVER.....	53
4.5 LAN DRIVER .....	54
4.6 USB 3.0 DRIVER .....	55
4.7 INTEL® MANAGEMENT ENGINE.....	56
4.8 KEYPAD AP .....	57
4.9 WIRELESS LAN DRIVER .....	59
4.10 BLUETOOTH DRIVER .....	59
4.11 RFID DRIVER (OPTIONAL).....	60
4.12 3-IN-1 COMBO READER DRIVER (OPTIONAL).....	63
4.12.1 SCR Driver.....	63
4.12.2 MSR Driver .....	67
4.12.3 Fingerprint Reader Driver.....	70
4.13 BARCODE READER DRIVER (OPTIONAL) .....	72
<b>5 BIOS SETUP .....</b>	<b>74</b>
5.1 INTRODUCTION.....	75
5.1.1 Starting Setup.....	75
5.1.2 Using Setup .....	76
5.1.3 Getting Help.....	76

## POC-W22A-H81 Medical Panel PC

5.1.4 BIOS Menu Bar.....	77
5.2 MAIN.....	77
5.3 ADVANCED.....	79
5.3.1 ACPI Settings.....	80
5.3.2 RTC Wake Settings.....	81
5.3.3 Trusted Computing.....	82
5.3.4 CPU Configuration.....	83
5.3.5 SATA Configuration.....	85
5.3.6 USB Configuration.....	86
5.3.7 F81866 Super IO Configuration.....	87
5.3.7.1 Serial Port n Configuration.....	88
5.3.7.1.1 Serial Port 1 Configuration.....	88
5.3.7.1.2 Serial Port 2 Configuration.....	89
5.3.8 iWDD H/W Monitor.....	90
5.3.8.1 Smart Fan Mode Configuration.....	91
5.3.9 Serial Port Console Redirection.....	93
5.3.10 IEI Feature.....	94
5.4 CHIPSET.....	95
5.4.1 PCH-IO Configuration.....	95
5.4.1.1 PCI Express Configuration.....	96
5.4.1.2 PCH Azalia Configuration.....	97
5.4.2 System Agent (SA) Configuration.....	98
5.4.2.1 Graphics Configuration.....	99
5.4.2.2 Memory Configuration.....	101
5.5 BOOT.....	101
5.6 SECURITY.....	103
5.7 SAVE & EXIT.....	104
5.8 SERVER MANAGEMENT.....	105
5.8.1 System Event Log.....	106
5.8.2 BMC Network Configuration.....	108
<b>6 SYSTEM MAINTENANCE.....</b>	<b>109</b>
6.1 SYSTEM MAINTENANCE INTRODUCTION.....	110
6.2 ANTI-STATIC PRECAUTIONS.....	110
6.3 TURN OFF THE POWER.....	111

6.4 SYSTEM FAN REPLACEMENT .....	111
6.5 SO-DIMM REPLACEMENT.....	112
6.6 REINSTALLING THE BACK COVER .....	114
6.7 TROUBLESHOOTING.....	114
<b>A REGULATORY COMPLIANCE .....</b>	<b>115</b>
<b>B SAFETY PRECAUTIONS .....</b>	<b>120</b>
B.1 SAFETY PRECAUTIONS.....	121
<i>B.1.1 General Safety Precautions .....</i>	<i>121</i>
<i>B.1.2 Anti-static Precautions .....</i>	<i>122</i>
<i>B.1.3 Product Disposal .....</i>	<i>123</i>
<i>B.1.4 Classification.....</i>	<i>124</i>
B.2 MAINTENANCE AND CLEANING PRECAUTIONS .....	124
<i>B.2.1 Maintenance and Cleaning.....</i>	<i>124</i>
<i>B.2.2 Cleaning Tools.....</i>	<i>125</i>
<b>C EMC TEST SUMMARY .....</b>	<b>126</b>
<b>D BIOS MENU OPTIONS.....</b>	<b>132</b>
<b>E WATCHDOG TIMER.....</b>	<b>135</b>
<b>F HAZARDOUS MATERIALS DISCLOSURE.....</b>	<b>138</b>



# List of Figures

---

Figure 1-1: POC-W22A-H81 Medical Panel PC .....	2
Figure 1-2: Front View .....	5
Figure 1-3: Backlit Touch Buttons .....	6
Figure 1-4: Side View .....	7
Figure 1-5: Bottom Panel .....	8
Figure 1-6: Rear View .....	8
Figure 1-7: Dimensions (mm) .....	13
Figure 3-1: Back Cover Retention Screws .....	23
Figure 3-2: Remove the Back Cover .....	23
Figure 3-3: HDD Bracket Retention Screws .....	24
Figure 3-4: HDD Bottom Bracket Removal .....	24
Figure 3-5: HDD Retention Screws .....	25
Figure 3-6: HDD Cover Installation .....	25
Figure 3-7: HDD Installation .....	26
Figure 3-8: Handset Holder Retention Screws .....	27
Figure 3-9: Handset Installation .....	28
Figure 3-10: Handset Driver Folder .....	28
Figure 3-11: Handset Driver Installation .....	29
Figure 3-12: Allow API Access .....	29
Figure 3-13: Manage Program Access to Skype .....	30
Figure 3-14: Handle Installation .....	32
Figure 3-15: Insert Barcode Reader Set .....	33
Figure 3-16: Install and Connect Barcode Reader Set.....	33
Figure 3-17: Barcode Reader Button .....	34
Figure 3-18: Reading Light .....	34
Figure 3-19: 3-in-1 Combo Reader .....	35
Figure 3-20: USB Ports on the Side Panel .....	35
Figure 3-21: Combo Reader Installation .....	36
Figure 3-22: Secure the Combo Reader .....	36
Figure 3-23: RFID Program Location .....	37
Figure 3-24: IRFR-100 Icon .....	37

Figure 3-25: IRFR Screen.....	38
Figure 3-26: IRFR – Find Tags.....	38
Figure 3-27: IRFR – UIDs.....	39
Figure 3-28: AT/ATX Switch Location.....	40
Figure 3-29: Aligning Tabs on the Bottom Panel .....	41
Figure 3-30: Cable Cover Installation .....	41
Figure 3-31: Cable Cover Removal .....	42
Figure 3-32: Wall-mounting Bracket.....	43
Figure 3-33: Chassis Support Screws.....	45
Figure 3-34: Secure the Panel PC .....	45
Figure 3-35: Arm Mounting Retention Screw Holes.....	46
Figure 3-36: Arm Mounting.....	47
Figure 3-37: Stand Mounting with EZ Stand (MEDP-EZS-R10).....	48
Figure 3-38: Powering On the System.....	49
Figure 3-39: Reset Button Location.....	49
Figure 4-1: Available Drivers .....	51
Figure 4-2: Intel® Chipset Device Software Install Wizard.....	52
Figure 4-3: Intel® Graphics Driver InstallShield Wizard.....	53
Figure 4-4: Realtek HD Audio Driver InstallShield Wizard .....	54
Figure 4-5: LAN Driver Install Wizard .....	55
Figure 4-6: USB 3.0 Driver Install Wizard.....	56
Figure 4-7: Intel® ME Components Install Wizard .....	57
Figure 4-8: Keypad Setup Wizard .....	58
Figure 4-9: Keypad AP .....	58
Figure 4-10: Wireless LAN InstallShield Wizard.....	59
Figure 4-11: Bluetooth Driver InstallShield Wizard.....	60
Figure 4-12: Device Manager - Update Driver Software.....	61
Figure 4-13: Update Driver Software Window .....	61
Figure 4-14: Browse for Driver Software Window .....	62
Figure 4-15: Driver Installation Complete Window .....	62
Figure 4-16: Device Manager Window–RFID Devices.....	63
Figure 4-17: Device Manager - Update Driver Software.....	64
Figure 4-18: Update Driver Software Window .....	64
Figure 4-19: Browse for Driver Software Window .....	65
Figure 4-20: Installing Driver Window .....	65

## POC-W22A-H81 Medical Panel PC

Figure 4-21: Driver Installation Complete Window .....	66
Figure 4-22: Device Manager Window–SCR Device.....	66
Figure 4-23: Device Manager - Update Driver Software.....	67
Figure 4-24: Update Driver Software Window .....	68
Figure 4-25: Browse for Driver Software Window .....	68
Figure 4-26: Installing Driver Window .....	69
Figure 4-27: Driver Installation Complete Window .....	69
Figure 4-28: Device Manager Window–MSR Device .....	70
Figure 4-29: Fingerprint Reader Driver Folder .....	71
Figure 4-30: Fingerprint Reader Driver InstallShield Wizard .....	71
Figure 4-31: Barcode Reader Driver Folder .....	72
Figure 4-32: Barcode Reader Driver Installation .....	72
Figure 4-33: Device Manager Window–Barcode Reader Device .....	73
Figure 6-1: System Fan Retention Screws and Connector .....	111
Figure 6-2: SO-DIMM Slot Cover Location .....	112
Figure 6-3: SO-DIMM Locations .....	113
Figure 6-4: SO-DIMM Installation .....	113

# List of Tables

---

Table 1-1: Model Variations .....	3
Table 1-2: Touch Button Functions .....	6
Table 1-3: System Specifications.....	12
Table 3-1: Handset Button Functions.....	31
Table 3-2: RS-232/422/485 Serial Port Pinouts .....	40
Table 5-1: BIOS Navigation Keys .....	76

# List of BIOS Menus

BIOS Menu 1: Main .....	78
BIOS Menu 2: Advanced .....	79
BIOS Menu 3: ACPI Settings .....	80
BIOS Menu 4: RTC Wake Settings .....	81
BIOS Menu 5: Trusted Computing .....	82
BIOS Menu 6: CPU Configuration .....	83
BIOS Menu 7: SATA Configuration .....	85
BIOS Menu 8: USB Configuration .....	86
BIOS Menu 9: F81866 Super IO Configuration .....	87
BIOS Menu 10: Serial Port n Configuration Menu .....	88
BIOS Menu 11: iWDD H/W Monitor .....	90
BIOS Menu 12: Smart Fan Mode Configuration .....	91
BIOS Menu 13: Serial Port Console Redirection .....	93
BIOS Menu 14: IEI Feature .....	94
BIOS Menu 15: Chipset .....	95
BIOS Menu 16: PCH-IO Configuration .....	95
BIOS Menu 17: PCI Express Configuration .....	96
BIOS Menu 18: PCH Azalia Configuration .....	97
BIOS Menu 19: System Agent (SA) Configuration .....	98
BIOS Menu 20: Graphics Configuration .....	99
BIOS Menu 21: Memory Configuration .....	101
BIOS Menu 22: Boot .....	101
BIOS Menu 23: Security .....	103
BIOS Menu 24: Save & Exit .....	104
BIOS Menu 25: Server Management .....	105
BIOS Menu 26: System Event Log .....	106
BIOS Menu 27: BMC Network Configuration .....	108





Chapter

1

# Introduction

---

## 1.1 Overview



**Figure 1-1: POC-W22A-H81 Medical Panel PC**

The POC-W22A-H81 is a 4<sup>th</sup> generation Intel® Core™/Pentium® CPU powered medical-grade panel PC with a rich variety of functions and peripherals. All POC-W22A-H81 models are designed for easy and simplified integration into point-of-care (POC) applications.

An Intel® Core™ i5/Core™ i3/Pentium® processor coupled with the Intel® H81 chipset delivers optimal memory, graphics, and peripheral I/O support. The system comes with 4.0 GB of pre-installed DDR3 SO-DIMMs and supports a maximum of 16.0 GB ensuring smooth data throughputs with reduced bottlenecks and fast system access. Dual display support is provided via the HDMI port or the VGA port.

One RS-232/422/485 serial port, two USB 3.0 ports and four USB 2.0 ports provide simplified connectivity to a variety of external peripheral devices. Wi-Fi 802.11a/b/g/n/ac high speed wireless and two RJ-45 GbE connectors allow for smooth connection of the system to an external LAN. The system also equips with a SATA interface, supporting both SATA HDD and SSD.

In addition, the POC-W22A-H81 features Intelligent Platform Management Interface 2.0 (IPMI 2.0) that helps lower the overall costs of server management by enabling users to maximize IT resource, save time and manage multiple systems. The POC-W22A-H81 supports IPMI 2.0 through the optional iRIS-2400 module.

## POC-W22A-H81 Medical Panel PC



### NOTE:

The POC-W22A-H81 medical panel PC is intended to be used to review and update electronic medical records (EMR) in hospital information system (HIS). The device is not suitable for diagnosis display.

Equipment connected to analog or digital interfaces must comply with the respective IEC Standards (e.g. IEC 60950 for data processing equipment and IEC 60601-1 for medical equipment). Furthermore all configurations shall comply with the current version of the standard for SYSTEMS IEC 60601-1-1. Everybody who connects additional equipment to the signal input part or signal output part configure a medical system, and is therefore responsible that the system complies with current version of the requirements of the system standard IEC 60601-1-1. If in doubt, consult the technical service department or your local representative.

## 1.2 Model Variations

There are three models in the POC-W22A-H81 series. All models are preinstalled with two 2 GB DDR3 memory modules and an 802.11a/b/n/ac Wi-Fi module. The model numbers and model variations are listed below.

Model	CPU
<b>POC-W22A-H81-P/PC/4G</b>	Intel® Pentium® G3320TE (dual-core, 2.3 GHz, TDP 35 W)
<b>POC-W22A-H81-i3/PC/4G</b>	Intel® Core™ i3-4330TE (dual-core, 2.4 GHz, TDP 35 W)
<b>POC-W22A-H81-i5/PC/4G</b>	Intel® Core™ i5-4570TE (dual-core, 2.7 GHz, TDP 35 W)

**Table 1-1: Model Variations**

## 1.3 Features

The POC-W22A-H81 features are listed below:


- 21.5" (16:9) flat-bezel LCD with LED backlight
- Anti-bacteria cover
- 4th generation Intel® Core™ i5, Core™ i3 or Pentium® processor
- Pre-installed with 4 GB of DDR3 SO-DIMM memory (system max. 16 GB)
- Projected capacitive type touchscreen
- Two PCIe GbE RJ-45 connectors and Wi-Fi 802.11a/b/g/n/ac high speed wireless
- Two internal 3 W speakers
- Four USB 2.0 ports and two USB 3.0 ports
- 12 V – 28 V wide range DC power input
- Optional internal 54 W battery module
- One RS-232/422/485 serial port by D-sub 9 connector
- IP 65 compliant front panel
- Optional Mifare RFID reader
- Optional 3-in-1 card reader (supports magnetic stripe card, smart card and fingerprint)
- Optional VoIP phone handset
- Optional handle module with 1D/2D barcode scanner and reading light



## POC-W22A-H81 Medical Panel PC

### 1.4 Front Panel

The front side of the POC-W22A-H81 is a flat-bezel panel with a TFT LCD screen surrounded by a PC/ABS plastic frame (**Figure 1-2**).



**Figure 1-2: Front View**

If the POC-W22A-H81 is installed with an optional UPS battery module, there will be a battery LED indicator located on the front panel. The UPS battery can provide backup power for 30 minutes in idle mode, and the maximum discharge voltage is 3.7 A. The POC-W22A-H81 will shutdown when the discharge voltage is over 3.7 A. The status descriptions of the battery LED indicator are listed below.

- **Solid amber:** the battery is being charged
- **Solid blue:** the system is using battery power or the battery is fully charged.
- **Blinking amber:** the battery is low (<30%). The battery indicator blinks regularly then turns off after 5 minutes.

Detail information about the UPS battery is described in **Section 3.6**.

**1.4.1 Backlit Touch Buttons**

The front panel of the POC-W22A-H81 contains several backlit touch buttons that control audio volume, LCD brightness and some other system components.



**Figure 1-3: Backlit Touch Buttons**

The following table describes the function of each button.

Button	Function
	LCD on/off (the touch buttons blink when the LCD is turned off)
	-: Brightness down or lock/unlock OSD (with function key) (minimum brightness: 5%) +: Brightness up (maximum brightness: 100%)
	-: Volume down +: Volume up
	Touch lock for cleaning: long-press for 3 seconds to lock or unlock the touch function of the screen. The touch buttons blink when the touch function is locked. The lock will be automatically disabled after 2 minutes.
	Lock or unlock the touch buttons: long-press for 3 seconds to lock or unlock the touch buttons. The touch buttons blink when the touch buttons are locked.
	Power on: long-press for 3 seconds. All touch buttons blink 3 times, and the system starts to boot. Power off: long-press for 3 seconds. All touch buttons blink 3 times, and the system starts to shut down.
Note: Press the touch button for at least one second to activate it.	

**Table 1-2: Touch Button Functions**

## POC-W22A-H81 Medical Panel PC

### 1.5 Side Panels

The left side panel has two USB 2.0 ports which are protected by a waterproof cover. The right side panel has a space for I/O interface expansion by installing a PCIe Mini I/O card.





Figure 1-4: Side View

### 1.6 Bottom Panel

The bottom panel of the POC-W22A-H81 has the following connectors and switches (Figure 1-5):

- 1 x 12 V ~ 28 V DC input power jack
- 1 x RS-232/422/485 DB-9 connector
- 1 x Barcode reader RJ-11connector
- 2 x GbE RJ-45 connector
- 2 x USB 3.0 connector
- 2 x USB 2.0 connector
- 1 x HDMI connector
- 1 x VGA connector
- 1 x AT/ATX switch
- 1 x Power button
- 1 x Reset button


**POC-W22A-H81 Medical Panel PC**



**Figure 1-5: Bottom Panel**

**1.7 Rear Panel**

The rear panel provides access to retention screw holes that support VESA mounting. See **Figure 1-6**.



**Figure 1-6: Rear View**



## POC-W22A-H81 Medical Panel PC

### 1.8 System Specifications

The technical specifications for the POC-W22A-H81 systems are listed in **Table 1-3**.


Specification	POC-W22A-H81
<b>CPU</b>	Intel® Pentium® G3320TE (dual-core, 2.3 GHz, TDP 35 W) Intel® Core™ i3-4330TE (dual-core, 2.4 GHz, TDP 35 W) Intel® Core™ i5-4570TE (dual-core, 2.7 GHz, TDP 35 W)
<b>Chipset</b>	Intel® H81
<b>Memory</b>	Two 204-pin 1600/1333 MHz dual-channel DDR3 SO-DIMM slots preinstalled with two 2 GB SDRAM (system max. 16 GB)
<b>Storage</b>	One PCIe Mini card slot for mSATA module installation (ATO) One 2.5" SATA HDD bay
<b>Auto-dimming</b>	Built-in ambient light sensor for panel brightness adjustment
<b>LCD and Touchscreen</b>	
<b>LCD Size</b>	21.5" (16:9)
<b>Max. Resolution</b>	1920 (W) x 1080 (H)
<b>Brightness (cd/m<sup>2</sup>)</b>	250
<b>Contrast Ratio</b>	1000:1
<b>LCD Color</b>	16.7M
<b>Pixel Pitch (mm)</b>	0.24825 (H) x 0.24825 (V)
<b>Viewing Angle (H-V)</b>	170° / 160°
<b>Backlight MTBF</b>	30,000 hrs
<b>Backlight</b>	LED
<b>Touchscreen</b>	Projected capacitive type with USB interface
<b>Surface Hardness</b>	6H



<b>Network Connection</b>	
<b>Wireless</b>	One pre-installed wireless LAN module (half-size PCIe Mini card) supports 802.11a/b/g/n/ac
<b>LAN</b>	Dual GbE connector
<b>Audio</b>	
<b>Audio</b>	Realtek ALC892 HD Audio codec
<b>Internal Speaker</b>	Two 3 W speakers
<b>Camera</b>	2-megapixel with auto focus and digital microphone
<b>Handset</b>	Optional VoIP phone handset
<b>Optional Features</b>	
<b>Battery</b>	Optional 3550 mA, 54 W internal Li-Ion battery with LED indicator
<b>Remote Management</b>	One iRIS-2400 slot supports iRIS remote management module
<b>RFID Reader</b>	Optional Mifare 13.56 MHz card reader (with LED indicator)
<b>Card Reader</b>	Optional 3-in-1 card reader supports magnetic stripe card, smart card and fingerprint
<b>Handle</b>	Optional handle with/without 1D/2D barcode reader and reading light
<b>Connectors</b>	
<b>I/O Ports</b>	<p>1 x 12 V ~ 28 V DC input jack</p> <p>1 x Barcode reader connector (RJ-11)</p> <p>1 x HDMI connector</p> <p>1 x RS-232/422/485 serial port (DB-9 connector)</p> <p>2 x GbE LAN (RJ-45 connector)</p> <p>2 x USB 3.0 connectors</p> <p>4 x USB 2.0 connectors</p>

## POC-W22A-H81 Medical Panel PC

	1 x VGA connector	
<b>Buttons, Switches and Indicators</b>		
<b>Backlit Touch Buttons</b>	Six touch buttons (LCD on/off, brightness up, brightness down, volume up, volume down, lock/unlock touch function)	
<b>Buttons &amp; Switches</b>	1 x Power button 1 x AT/ATX switch 1 x Reset button	
<b>LED Indicators</b>	2-light battery status LED indicator (optional) RFID LED indicator (optional)	
<b>Physical</b>		
<b>Construction Material</b>	PC+ABS plastic with anti-bacterial material	
<b>VESA Mount</b>	100 mm x 100 mm	
<b>Dimensions (W x H x D)</b>	542.5 mm x 349.5 mm x 52 mm	
<b>Net Weight</b>	7.3 kg	
<b>Environment</b>		
<b>Transportation/Storage</b>	<b>Temperature</b>	-20°C ~ 60°C
	<b>Humidity</b>	10% ~ 95% (non-condensing)
	<b>Pressure</b>	700 hPa ~ 1060 hPa
<b>Operating</b>	<b>Temperature</b>	0°C ~ 40°C
	<b>Humidity</b>	10% ~ 95% (non-condensing)
	<b>Pressure</b>	700 hPa ~ 1060 hPa
<b>Vibration</b>	1G	
<b>Shock</b>	Operating Shock: 5G peak acceleration (11ms duration) Non-Operating Shock: 15G peak acceleration (11ms duration)	


<b>IP Level</b>	IP 65 compliant front panel
<b>Safety/EMC</b>	CE, FCC class B part 18, UL 60601-1, EN 60601-1
<b>Power</b>	
<b>Power Supply</b>	120 W medical-grade power adapter (FSP PMP120-13-2)
	Input: 100 V AC ~ 240 V AC, 47 Hz ~ 63 Hz, 1.4 A ~ 0.6 A
	Output: 120 W Max., 19 V  6.32 A
<b>Power Requirement</b>	12 V ~ 28 V DC

**Table 1-3: System Specifications**

**POC-W22A-H81 Medical Panel PC**

**1.9 Dimensions**

The POC-W22A-H81 dimensions are shown below.



**Figure 1-7: Dimensions (mm)**

Chapter

**2**

# Unpacking

---



## POC-W22A-H81 Medical Panel PC

### 2.1 Unpacking

To unpack the medical panel PC, follow the steps below:



#### **WARNING!**

The front side LCD screen has a protective plastic cover stuck to the screen. Only remove the plastic cover after the medical panel PC has been properly installed. This ensures the screen is protected during the installation process.

---

- Step 1:** Use box cutters, a knife or a sharp pair of scissors that seals the top side of the external (second) box.
- Step 2:** Open the external (second) box.
- Step 3:** Use box cutters, a knife or a sharp pair of scissors that seals the top side of the internal (first) box.
- Step 4:** Lift the monitor out of the boxes.
- Step 5:** Remove both polystyrene ends, one from each side.
- Step 6:** Pull the plastic cover off the medical panel PC.
- Step 7:** Make sure all the components listed in the packing list are present.






## 2.2 Packing List





**NOTE:**

If any of the components listed in the checklist below are missing, do not proceed with the installation. Contact the IEI reseller or vendor the POC-W22A-H81 was purchased from or contact an IEI sales representative directly by sending an email to [sales@ieiworld.com](mailto:sales@ieiworld.com).

The POC-W22A-H81 medical panel PC is shipped with the following components:




Quantity	Item	Image
1	POC-W22A-H81 medical panel PC	
1	Medical-grade power adapter (120 W, 19 V DC output) (P/N: 63040-010120-010-RS)	
1	Power cord (P/N: 32702-000200-100-RS)	
4	Pan-head screw (M3*5) for HDD installation (P/N: 44043-030051-RS)	
1	Quick Installation Guide	

**POC-W22A-H81 Medical Panel PC**

1	Utility CD	
1	One Key Recovery CD	





**2.3 Optional Items**

The following are optional components which may be separately purchased:

Item and Part Number	Image
3-in-1 reader (smart card / magnetic card / fingerprint) (P/N: MEDP-CR-R10)	
Handle (P/N: MEDP-HD-R10)	
Handle with 1D/2D barcode reader and reading light (P/N: MEDP-HD-BR-R10)	

Item and Part Number	Image
Handset (USB interface) and holder (P/N: MEDP-HS-R10)	 A white, vertical handset with a coiled cable and a mounting bracket.
Cable cover (P/N: POCP-CC01-R10)	 A long, thin, white plastic cable cover.
EZ Stand with cabling hole, VESA 100 (P/N: MEDP-EZS-R10)	 A white, adjustable stand with a VESA 100 mounting hole and a cabling hole.
VESA 100 wall mount kit (four M3*6 screws included) (P/N: AFLWK-19B)	 A silver metal VESA 100 wall mount kit.
Arm (P/N: ARM-31-RS)	 A yellow, adjustable arm for mounting a monitor.
Stand (P/N: STAND-A21-R10)	 A black, adjustable stand for a monitor.

## POC-W22A-H81 Medical Panel PC

Item and Part Number	Image
iRIS module (assemble-to-order) (P/N: iRIS-2400-R10)	
Trusted platform module (assemble-to-order) (P/N: MEDP-TPM-R10)	
14.8 V 3550 mAH Li-ion battery (Getac BAT003-4S1P3550-0, assemble-to-order) (P/N: MEDP-BAT-R10)	
Mifare RFID reader compliant with ISO 14443A, ISO 14443B and ISO 15693 protocols (assemble-to-order) (P/N: MEDP-MF-RFID-R10)	

If any of these items are missing or damaged, contact the distributor or sales representative immediately.



Chapter

**3**

# Installation

---

### 3.1 Anti-static Precautions

---

**WARNING:**

Failure to take ESD precautions during the maintenance of the POC-W22A-H81 may result in permanent damage to the POC-W22A-H81 and severe injury to the user.

---

Electrostatic discharge (ESD) can cause serious damage to electronic components, including the POC-W22A-H81. Dry climates are especially susceptible to ESD. It is therefore critical that whenever the POC-W22A-H81 is accessed internally, or any other electrical component is handled, the following anti-static precautions are strictly adhered to.

- ***Wear an anti-static wristband:*** Wearing a simple anti-static wristband can help to prevent ESD from damaging the board.
- ***Self-grounding:*** Before handling the board, touch any grounded conducting material. During the time the board is handled, frequently touch any conducting materials that are connected to the ground.
- ***Use an anti-static pad:*** When configuring the POC-W22A-H81, place it on an anti-static pad. This reduces the possibility of ESD damaging the POC-W22A-H81.
- ***Only handle the edges of the PCB:*** When handling the PCB, hold the PCB by the edges.

### 3.2 Installation Precautions

When installing the medical panel PC, please follow the precautions listed below:

- **Manufacturer authorization:** Do not modify this equipment without authorization of manufacturer.
- **Certified Engineers:** Only certified engineers should install and modify the hardware settings.

- **Power turned off:** When installing the medical panel PC, make sure the power is off. Failing to turn off the power may cause severe injury to the body and/or damage to the system.
- **Anti-static Discharge:** If a user open the rear panel of the medical panel PC, to configure the jumpers or plug in added peripheral devices, ground themselves first and wear an anti-static wristband.

### 3.3 Installation and Configuration Steps

The following installation steps must be followed.

**Step 1:** Unpack the medical panel PC.

**Step 2:** Install a HDD.

**Step 3:** Install necessary accessories (handset, handle or combo reader)

**Step 4:** Configure the system.

**Step 5:** Connect peripheral devices to the medical panel PC.

**Step 6:** Mount the medical panel PC.

### 3.4 Removing the Back Cover

To access the POC-W22A-H81 internally the back cover must be removed. To remove the back cover, please follow the steps below.

**Step 1:** Remove the two retention screws from the back cover (**Figure 3-1**).

## POC-W22A-H81 Medical Panel PC



**Figure 3-1: Back Cover Retention Screws**

**Step 2:** The back cover is attached to the chassis through a magnetic force. To remove the back cover, grasp the back cover finger grips (as shown in **Figure 3-2**) and lift the cover up.




**Figure 3-2: Remove the Back Cover**

### 3.5 HDD Installation

To install the HDD into the system, please follow the steps below:


**Step 1:** Remove the back cover. See **Section 3.4** above.

**Step 2:** Remove the two HDD bracket retention screws (**Figure 3-3**) and lift the HDD bracket off the panel PC.



**Figure 3-3: HDD Bracket Retention Screws**

**Step 3:** Slide the bottom bracket to remove it from the HDD cover (**Figure 3-4**).



**Figure 3-4: HDD Bottom Bracket Removal**




## POC-W22A-H81 Medical Panel PC

**Step 4:** Attach the bottom bracket to the HDD. To do this, align the four retention screw holes in the bottom of the HDD bracket with the retention screw holes on the bottom of the HDD. Insert four retention screws (M3\*5) into the bracket (**Figure 3-5**).



**Figure 3-5: HDD Retention Screws**

**Step 5:** Slide the bottom bracket into the HDD cover (**Figure 3-6**).



**Figure 3-6: HDD Cover Installation**

**Step 6:** Install the HDD into the POC-W22A-H81 by inserting and connecting the HDD to the SATA connectors. Then, secure the HDD bracket to the chassis with the previously removed retention screws. See **Figure 3-7**.




Figure 3-7: HDD Installation

**Step 7:** Replace the back cover and secure it using two retention screws.

### 3.6 UPS Battery

An optional UPS battery can be assembled to order to provide backup power for up to 30 minutes in case of a power interruption. The battery LED indicator is located on the front panel. The LED status descriptions are described in **Section 1.4**.



#### CAUTION:

The following precautions should be followed when using the battery:

- Do not use battery power to boot up the POC-W22A-H81.
  - Charge the battery with a voltage higher than 16.8 V.
  - Charging time (from 30% to 100%):
    - System off: 3 hours
    - System on: 2 hours
  - Risk of explosion if battery is replaced by an incorrect type. Only certified engineers should replace the on-board battery.
  - Dispose of used batteries according to instructions and local regulations.
-

## POC-W22A-H81 Medical Panel PC

### 3.7 Handset Installation (Optional)

An optional phone handset can be installed on the side of the POC-W22A-H81 to make VoIP calls. To install the handset and the holder, please follow the instruction below.

- Step 1:** Locate the two retention screw holes for installing the handset holder on the rear panel.
- Step 2:** Secure the handset holder with the POC-W22A-H81 by two retention screws (M3\*8, flat head).



**Figure 3-8: Handset Holder Retention Screws**

- Step 3:** Plug the handset cord into a USB connector on the bottom panel.
- Step 4:** Place the handset in the holder.




**Figure 3-9: Handset Installation**

### 3.7.1 Using VoIP Handset

The VoIP handset is designed for Skype. To use the handset to place or receive a call via Skype, please follow the steps below.

**Step 1:** Install the Skype program (<http://www.skype.com/en/>).

**Step 2:** Select **Other** from the list of the driver CD. Double click the setup file in the **POCP-W22A-HS-R10\_U-2000H** folder to install the handset driver (Figure 3-10).



**Figure 3-10: Handset Driver Folder**



## POC-W22A-H81 Medical Panel PC

**Step 3:** Follow the step-by-step instruction of the installation wizard (Figure 3-11) to install the handset driver.




Figure 3-11: Handset Driver Installation

**Step 4:** Launch Skype. Press the **Allow access** button on Skype (Figure 3-12) to allow handset API access.





Figure 3-12: Allow API Access

API access can also be managed through **Tools** → **Options** → **Advanced settings** in Skype. See Figure 3-13.









**Figure 3-13: Manage Program Access to Skype**

**Step 5:** The user can now use Skype via the handset. The function description of each button on the handset is listed in the following table.

## POC-W22A-H81 Medical Panel PC

Button	Function
	<b>LED indicator:</b> Clear - no call activity Blinking green - incoming call ringing Steady green - active call
	<b>Hot key:</b> No call activity - launch Skype and select menus * long-press the hot key for 2 seconds to turn off Skype Ringing - terminate the incoming call Active call - mute or unmute the handset microphone
	No call activity - scroll up through incoming and outgoing call history Active call - handset speaker volume up
	No call activity - scroll down through incoming and outgoing call history Active call - handset speaker volume down
	Place, answer or hang up a call.

**Table 3-1: Handset Button Functions**

### 3.8 Handle Installation (Optional)

An optional handle can be installed on the POC-W22A-H81 for the user to easily adjust the viewing angle and the position of the POC-W22A-H81. To install the handle, please follow the instruction below.

**Step 1:** Locate the retention screw holes for installing the handle on the rear panel. If a cable cover is installed on the rear panel, please remove it first.

**Step 2:** Secure the handle with the POC-W22A-H81 by inserting eight retention screws (M3\*6L, flat head).



**Figure 3-14: Handle Installation**

### **3.8.1 Barcode Reader Installation**

The optional handle may come with a barcode reader set which also contains a reading light with three levels of brightness. To install the barcode reader set, please follow the instruction below.

- Step 1:** Follow the instruction described above to install the handle.
- Step 2:** Insert the barcode reader set into the slot in the center of the handle. To be able to insert the barcode reader, the side with barcode reader must face toward the right side of the POC-W22A-H81 as shown in **Figure 3-15**.

## POC-W22A-H81 Medical Panel PC



**Figure 3-15: Insert Barcode Reader Set**

**Step 3:** Push the barcode reader set all the way down, and then rotate the barcode reader anti-clockwise to a proper position (**Figure 3-16**).

**Step 4:** Connect the barcode reader cable to the RJ-11 connector on the bottom panel of the POC-W22A-H81 (**Figure 3-16**).



**Figure 3-16: Install and Connect Barcode Reader Set**

**Step 5:** Install the driver for the barcode reader by following the instructions described in **Section 4.13**.

**Step 6:** After driver installation is complete, push the barcode reader button to trigger the barcode reader.





**Figure 3-17: Barcode Reader Button**



**WARNING:**

Do not stare into beam of the laser light. The human eye can be damaged. Avoid that the laser beam hits reflective surfaces such as mirrors, etc. Any changes at the device are forbidden these could cause a dangerous laser light.

### 3.8.2 Reading Light

The barcode reader is also equipped with a reading light with 3-level of brightness. Push the reading light button to turn on or to toggle illumination brightness levels.



**Figure 3-18: Reading Light**



## POC-W22A-H81 Medical Panel PC

### 3.9 3-in-1 Combo Reader Installation (Optional)

The 3-in-1 combo reader is an optional item for the POC-W22A-H81. The combo reader combines fingerprint reader, smart card reader (SCR) and magnetic stripe reader (MSR) into one compact device.




Figure 3-19: 3-in-1 Combo Reader

To install the combo reader to the POC-W22A-H81, please follow the steps below.

**Step 1:** Open the USB port cover on the side panel of the POC-W22A-H81. Secure the cover by rotating the cover and inserting the stub into the hole on the side panel.



Figure 3-20: USB Ports on the Side Panel

**Step 2:** Align the USB connectors on the reader with the two USB connectors on the side panel of the POC-W22A-H81.

**Step 3:** Insert and connect the USB connectors to install the combo reader.

See **Figure 3-21**.



**Figure 3-21: Combo Reader Installation**

**Step 4:** Secure the combo reader to the system by inserting two retention screws (M3x5L) into the rear panel of the POC-W22A-H81 and tightening them.



**Figure 3-22: Secure the Combo Reader**

## POC-W22A-H81 Medical Panel PC

**Step 5:** Install the drivers for these three readers by following the instructions described in **Chapter 4**:

- Section 4.12.1: SCR Driver
- Section 4.12.2: MSR Driver
- Section 4.12.3: Fingerprint Reader Driver


### 3.10 Using RFID Reader (Optional)

The POC-W22A-H81 may come with an optional RFID reader pre-installed inside the bottom of the front panel. To use the RFID reader, follow the steps below.

**Step 1:** Install the RFID driver (refer to **Section 4.11**).

**Step 2:** Locate the **IRFR-100.exe** file in the following folder of the driver CD:

\\Docs\10.Other\POCP-MF-RFID-R10\RFID\D490. Copy the **IRFR-100.exe** program to the desktop.




**Figure 3-23: RFID Program Location**

**Step 3:** Double click the **IRFR-100** icon on the desktop.




**Figure 3-24: IRFR-100 Icon**

**Step 4:** The IRFR-100 window appears (Figure 3-25).




**Figure 3-25: IRFR Screen**


**Step 5:** Select the **Find tags** tab and click the **Run** button to enable the RFID reader (Figure 3-26).



**Figure 3-26: IRFR – Find Tags**

## POC-W22A-H81 Medical Panel PC

**Step 6:** Place an RFID card near the RFID reader  on the bottom of the front panel (**Figure 1-2**) then remove it. The card number will be shown in the UIDs column (**Figure 3-27**).



**Figure 3-27:** IRFR – UIDs



### **NOTE:**

Please refer to the IRFR-100 user guide in the driver CD (IRFR-100\_AP\_UserGuide.pdf) for detailed instruction on how to use the IRFR-100.



### 3.11 RS-232/422/485 Serial Port Connection

The bottom panel of the POC-W22A-H81 has one D-sub 9 male connectors for RS-232/422/485 connection. The serial communication mode selection can be made through the BIOS options. Please refer to **Section 5.3.7.1.1** for detail information. The pinouts of the D-sub 9 connector are listed below.

Pin	RS-232	RS-422	RS-485
1	DCD	TXD422-	TXD485-
2	RX	TXD422+	TXD485+
3	TX	RXD422+	
4	DTR	RXD422-	
5	GND		
6	DSR		
7	RTS		
8	CTS		
9	RI		



**Table 3-2: RS-232/422/485 Serial Port Pinouts**

### 3.12 AT/ATX Mode Selection

AT or ATX power mode can be used on the POC-W22A-H81. The selection is made through an AT/ATX switch located on the bottom panel (**Figure 3-28**).



**Figure 3-28: AT/ATX Switch Location**

#### 3.12.1 AT Power Mode

With the AT mode selected, the power is controlled by a central power unit rather than a power switch. The POC-W22A-H81 panel PC turns on automatically when the power is connected.

## POC-W22A-H81 Medical Panel PC


### 3.12.2 ATX Power Mode

With the ATX mode selected, the POC-W22A-H81 panel PC goes in a standby mode when it is turned off. The panel PC can be easily turned on via network or a power switch in standby mode.

### 3.13 Cable Cover Installation


An optional cable cover can be installed on the POC-W22A-H81 for the user to easily manage cables. To install the cable cover, please follow the instruction below.

**Step 1:** Align the three tabs on the rear of the cable cover with the slots on the bottom panel of the POC-W22A-H81 (**Figure 3-29**). Then, insert the tabs into the slots.



**Figure 3-29: Aligning Tabs on the Bottom Panel**

**Step 2:** Push the two tabs on the cable cover (as shown in **Figure 3-30**) into the slots on the rear panel of the POC-W22A-H81 one by one. More strength is required to push the tabs into the slots.



**Figure 3-30: Cable Cover Installation**

**Step 3:** To remove the cable cover, push the two tabs inwards to release the cover (Figure 3-31), and lift the cover from the POC-W22A-H81.



**Figure 3-31: Cable Cover Removal**

### 3.14 Mounting the System

The methods of mounting the POC-W22A-H81 are listed below.

- Wall mounting
- Arm mounting
- Stand mounting

The mounting methods are described below.



**WARNING:**

Use suitable mounting apparatus and be sure to secure the screws of the mounting apparatus tightly to avoid risk of injury.


---

## POC-W22A-H81 Medical Panel PC

### 3.14.1 Wall Mounting

To mount the medical panel PC onto the wall, please follow the steps below.

- Step 1:** Select the location on the wall for the wall-mounting bracket.
- Step 2:** Carefully mark the locations of the four screw holes in the bracket on the wall.
- Step 3:** Drill four pilot holes at the marked locations on the wall for the bracket retention screws.
- Step 4:** Align the wall-mounting bracket screw holes with the pilot holes.
- Step 5:** Secure the mounting-bracket to the wall by inserting the retention screws into the four pilot holes and tightening them (**Figure 3-32**).




**Figure 3-32: Wall-mounting Bracket**

**Step 6:** Insert the four monitor mounting screws provided in the wall mount kit into the four screw holes on the rear panel of the medical panel PC and tighten until the screw shank is secured against the rear panel (**Figure 3-33**).

**WARNING:**

Please use the M4 screws provided in the wall mount kit for the rear panel. If the screw is missing, the thread depth of the replacement screw should be not more than 4 mm.



---

**Step 7:** Align the mounting screws on the monitor rear panel with the mounting holes on the bracket.

**Step 8:** Carefully insert the screws through the holes and gently pull the monitor downwards until the monitor rests securely in the slotted holes (**Figure 3-33**). Ensure that all four of the mounting screws fit snugly into their respective slotted holes.

**NOTE:**

In the diagram below the bracket is already installed on the wall.

---



POC-W22A-H81 Medical Panel PC




Figure 3-33: Chassis Support Screws

**Step 9:** Secure the panel PC by fastening the retention screw of the wall-mounting bracket (Figure 3-34).




Figure 3-34: Secure the Panel PC

### 3.14.2 Arm Mounting

The POC-W22A-H81 is VESA (Video Electronics Standards Association) compliant and can be mounted on an arm with a 100 mm interface pad. To mount the POC-W22A-H81 on an arm, please follow the steps below.

**Step 1:** The arm is a separately purchased item. Please correctly mount the arm onto the surface it uses as a base. To do this, refer to the installation documentation that came with the mounting arm.




#### NOTE:

When purchasing the arm please ensure that it is VESA compliant and that the arm has a 100 mm interface pad. If the mounting arm is not VESA compliant it cannot be used to support the POC-W22A-H81 medical panel PC.

---

**Step 2:** Once the mounting arm has been firmly attached to the surface, lift the panel PC onto the interface pad of the mounting arm.


**Step 3:** Align the retention screw holes on the mounting arm interface with those in the panel PC (**Figure 3-35**).



**Figure 3-35: Arm Mounting Retention Screw Holes**

## POC-W22A-H81 Medical Panel PC

- Step 4:** Secure the POC-W22A-H81 to the interface pad by inserting four retention screws through the mounting arm interface pad and into the POC-W22A-H81.




**Figure 3-36: Arm Mounting**

### 3.14.3 Stand Mounting

To mount the POC-W22A-H81 using the EZ stand mounting kit, please follow the steps below.

- Step 1:** Locate the VESA mount screw holes on the rear of the POC-W22A-H81 (Figure 3-35). This is where the bracket will be attached.
- Step 2:** Align the bracket with the screw holes.
- Step 3:** To secure the bracket to the POC-W22A-H81 insert the retention screws into the screw holes and tighten them.



**Figure 3-37: Stand Mounting with EZ Stand (MEDP-EZS-R10)**



**NOTE:**

If the EZ stand is mounted, the handle (MEDP-HD-R10 or MEDP-HD-BR-R10) can not be installed.

---

### 3.15 Powering On the System

---



**WARNING:**

To avoid risk of electric shock, this equipment must only be connected to supply mains with protective earth.

---



## POC-W22A-H81 Medical Panel PC



### CAUTION:

The FSP PMP120-13-2 power adapter came with the POC-W22A-H81 is a forming part of the medical device.

To power on the system, follow the steps below:

- Step 1:** Connect the power cord to the power adapter. Connect the other end of the power cord to a power source.
- Step 2:** Connect the power adapter to the power connector of the POC-W22A-H81.
- Step 3:** Locate the power button on the I/O panel.
- Step 4:** Short press the power button to turn on the POC-W22A-H81.



The user can also long-press the touch buttons  +  on the front panel for three seconds to power on the system (please refer to **Table 1-2**).



Figure 3-38: Powering On the System

### 3.16 Reset the System

The reset button enables user to reboot the system when the system is turned on. The reset button location is shown in **Figure 3-39**. Press the reset button to reboot the system.



Figure 3-39: Reset Button Location



Chapter

**4**

# Driver Installation

---

## POC-W22A-H81 Medical Panel PC

## 4.1 Available Software Drivers

**NOTE:**

The contents of the driver folder may vary throughout the life cycle of the product and is subject to change without prior notice. Visit the IEI website or contact technical support for the latest updates.

All the drivers for the POC-W22A-H81 are on the utility CD that came with the system. To install the drivers, please follow the steps below.

**Step 1:** Insert the CD into a CD drive connected to the system.

**Step 2:** The driver main menu with a list of available drivers appears (**Figure 4-1**).



**Figure 4-1: Available Drivers**


**Step 3:** Install all of the necessary drivers in this menu.

## 4.2 Intel® Chipset Driver

To install the chipset driver, please follow the steps below.

**Step 1:** Select **Chipset** from the list of the driver CD.

**Step 2:** Double click the setup file in the folder. The Intel® Chipset Device Software install wizard appears



**Figure 4-2: Intel® Chipset Device Software Install Wizard**

**Step 3:** Follow the step-by-step instruction of the installation wizard to install the driver.

## 4.3 Intel® Graphics Driver

To install the graphics driver, please follow the steps below.

**Step 1:** Select **VGA** from the list of the driver CD. Locate the driver setup file for the corresponding operating system.

**Step 2:** Double click the setup file in the folder. The **Intel® Graphics Driver InstallShield Wizard** appears (**Figure 4-3**).

## POC-W22A-H81 Medical Panel PC




Figure 4-3: Intel® Graphics Driver InstallShield Wizard

**Step 3:** Follow the step-by-step instruction of the installation wizard to install the graphics driver.

## 4.4 Audio Driver

To install the driver for the speaker and the microphone, please follow the steps below.

**Step 1:** Select **Audio** from the list of the driver CD.

**Step 2:** Double click the setup file in the folder. The **InstallShield Wizard** screen appears (**Figure 4-4**).



**Figure 4-4: Realtek HD Audio Driver InstallShield Wizard**

**Step 3:** Follow the step-by-step instruction of the installation wizard to install the HD Audio driver.

## 4.5 LAN Driver

To install the LAN driver, please follow the steps below.

**Step 1:** Select **LAN** from the list of the driver CD. Locate the driver setup file for the corresponding operating system.

- PROWin7\_Win8.1\_32\_18.7.exe: Windows 7 or 8.1 (32-bit)
- PROWin7\_Win8.1\_64\_18.7.exe: Windows 7 or 8.1 (64-bit)

**Step 2:** Double click the setup file in the folder. The **Install Wizard** screen appears (Figure 4-4).



## POC-W22A-H81 Medical Panel PC




Figure 4-5: LAN Driver Install Wizard


**Step 3:** Follow the step-by-step instruction of the installation wizard to install the Intel® Network Connection driver.

## 4.6 USB 3.0 Driver

To install the USB 3.0 driver, please follow the steps below.

**Step 1:** Select **USB 3.0** from the list of the driver CD. Locate the driver setup file.

**Step 2:** Double click the setup file. The **Intel® Installation Framework** window appears (Figure 4-3).



**Figure 4-6: USB 3.0 Driver Install Wizard**

**Step 3:** Follow the step-by-step instruction of the installation wizard to install the USB 3.0 driver.


## **4.7 Intel® Management Engine**

To install the Intel® Management Engine Components, please follow the steps below.

**Step 1:** Select **MEI** from the list of the driver CD. Locate the driver setup file.

**Step 2:** Double click the setup file. The **Intel® Installation Framework** window appears (Figure 4-3).

## POC-W22A-H81 Medical Panel PC



**Figure 4-7: Intel® ME Components Install Wizard**


**Step 3:** Follow the step-by-step instruction of the installation wizard to install the Intel® Management Engine Components.

## 4.8 Keypad AP

The Keypad AP is an OSD control tool developed by IEI. To install the Keypad AP, please follow the steps below.

**Step 1:** Select **Keypad AP** from the list of the driver CD.

**Step 2:** Double click the **KeypadAP H779 v1.0\_1.exe** file. The following window appears (**Figure 4-3**).



**Figure 4-8: Keypad Setup Wizard**

**Step 3:** Follow the step-by-step instruction of the installation wizard to install the Keypad AP.

**Step 4:** After the installation, the Keypad AP can be accessed by pressing the brightness up/down buttons or the volume up/down buttons on the bottom frame of the monitor. It allows users to control screen brightness, audio volume and auto-dimming function.



**Figure 4-9: Keypad AP**

## POC-W22A-H81 Medical Panel PC

### 4.9 Wireless LAN Driver


To install the wireless LAN driver, please follow the steps below.

**Step 1:** Select **WiFi** from the list of the driver CD. Locate the setup file in this folder:

**AZ\_RTL8821AE\_Win7\_Win8.X\_\_2012.9.0212.2014.**

**Step 2:** Double click the setup file in the folder. The InstallShield Wizard screen appears

(Figure 4-10).



**Figure 4-10: Wireless LAN InstallShield Wizard**

**Step 3:** Follow the step-by-step instruction of the installation wizard to install the

Wireless LAN driver.

### 4.10 Bluetooth Driver


To install the Bluetooth driver, please follow the steps below.

**Step 1:** Select **WiFi** from the list of the driver CD. Locate the setup file in this folder:

**BT\_805.806.806.0221.2014.**

**Step 2:** Double click the setup file in the folder. The InstallShield Wizard screen appears.





**Figure 4-11: Bluetooth Driver InstallShield Wizard**

**Step 3:** Follow the step-by-step instruction of the installation wizard to install the Bluetooth driver.

### **4.11 RFID Driver (Optional)**

To install the RFID driver, please follow the steps below.

**Step 1:** Open the Device Manager window. Long press or right click **USB <-> Serial**. Select **Update Driver Software** from the pop-up window.

# POC-W22A-H81 Medical Panel PC




Figure 4-12: Device Manager - Update Driver Software

**Step 2:** The **Update Driver Software** window appears. Select **Browse my computer for driver software** for driver software.





Figure 4-13: Update Driver Software Window


**Step 3:** The following window appears. Press/Click the **Browse** button to specify the RFID driver directory (\Docs\10.Other\POCP-MF-RFID-R10\RFID\D490). Then, press/click the **Next** button.



**Figure 4-14: Browse for Driver Software Window**

**Step 4:** The system starts installing the RFID driver.

**Step 5:** After the driver installation process is complete, a confirmation screen appears. Click **Close** to exit the program.



**Figure 4-15: Driver Installation Complete Window**

## POC-W22A-H81 Medical Panel PC

**Step 6:** Repeat **Step 1 ~ Step 5** to install the RFID driver again.

**Step 7:** The **Device Manager Window** now shows the installed RFID devices.




Figure 4-16: Device Manager Window – RFID Devices

### 4.12 3-in-1 Combo Reader Driver (Optional)

The drivers for the optional 3-in-1 combo reader are all located in the following folder of the driver CD: `\Docs\10.Other\POCP-W22A-CR-R10`. Please follow the instructions below to install the drivers.

#### 4.12.1 SCR Driver

Follow the steps below to install the SCR driver.

**Step 1:** Open the Device Manager window. Long press or right click **Singular VCOM Card Reader**. Select **Update Driver Software** from the pop-up window.




Figure 4-17: Device Manager - Update Driver Software

**Step 2:** The **Update Driver Software** window appears. Select **Browse my computer** for driver software.





Figure 4-18: Update Driver Software Window




## POC-W22A-H81 Medical Panel PC

**Step 3:** The following window appears. Press/Click the **Browse** button to specify the SCR driver directory (\Docs\10.Other\POCP-W22A-CR-R10\SCR). Then, press/click the **Next** button.




**Figure 4-19: Browse for Driver Software Window**

**Step 4:** The following window (**Figure 4-20**) appears as the driver is installed.




**Figure 4-20: Installing Driver Window**

**Step 5:** After the driver installation process is complete, a confirmation screen appears.  
Click **Close** to exit the program.



**Figure 4-21: Driver Installation Complete Window**

**Step 6:** The **Device Manager Window** now shows the installed SCR device.



**Figure 4-22: Device Manager Window – SCR Device**

## POC-W22A-H81 Medical Panel PC

### 4.12.2 MSR Driver

Follow the steps below to install the MSR driver.

**Step 1:** Open the Device Manager window. Long press or right click **Singular VCOM Card Reader**. Select **Update Driver Software** from the pop-up window.






Figure 4-23: Device Manager - Update Driver Software

**Step 2:** The **Update Driver Software** window appears. Select **Browse my computer for driver software**.



**Figure 4-24: Update Driver Software Window**

**Step 3:** The following window appears. Press/Click the **Browse** button to specify the MSR driver directory (\Docs\10.Other\POCP-W22A-CR-R10\MSR). Then, press/click the **Next** button.



**Figure 4-25: Browse for Driver Software Window**

**Step 4:** The following window (**Figure 4-26**) appears as the driver is installed.

## POC-W22A-H81 Medical Panel PC




Figure 4-26: Installing Driver Window

**Step 5:** After the driver installation process is complete, a confirmation screen appears.  
Click **Close** to exit the program.




Figure 4-27: Driver Installation Complete Window

**Step 6:** The **Device Manager Window** now shows the installed MSR device.






Figure 4-28: Device Manager Window – MSR Device

### 4.12.3 Fingerprint Reader Driver

Follow the steps below to install the fingerprint reader driver.

**Step 1:** Select **Other** from the list of the driver CD.

**Step 2:** The fingerprint reader driver is located in the following folder (**Figure 4-29**):  
**\\Docs\10.Other\POCP-W22A-CR-R10\Finger Printer**. Double click the **setup.exe** file in this folder to install the driver.

## POC-W22A-H81 Medical Panel PC




Figure 4-29: Fingerprint Reader Driver Folder

**Step 3:** The Egis ES603 WDM Driver welcome window appears.





Figure 4-30: Fingerprint Reader Driver InstallShield Wizard

**Step 4:** Follow the step-by-step instruction of the installation wizard to install the fingerprint reader driver.

### 4.13 Barcode Reader Driver (Optional)


To install the barcode reader driver, please follow the steps below.

**Step 1:** Select **Other** from the list of the driver CD. Double click the **Install\_x86.bat** file (or **Install\_x64.bat** for 64-bit OS) in the **POCP-W22A-HD-BR-R10** folder shown in **Figure 4-31** to install the barcode reader driver.



**Figure 4-31: Barcode Reader Driver Folder**

**Step 2:** The following window shows and starts installing the barcode reader driver. When the installation is complete, the window will close automatically.



**Figure 4-32: Barcode Reader Driver Installation**

**Step 3:** The **Device Manager Window** now shows the installed barcode reader device.

# POC-W22A-H81 Medical Panel PC




Figure 4-33: Device Manager Window – Barcode Reader Device

Chapter

**5**

# BIOS Setup

---



## 5.1 Introduction

A licensed copy of the BIOS is preprogrammed into the ROM BIOS. The BIOS setup program allows users to modify the basic system configuration. This chapter describes how to access the BIOS setup program and the configuration options that may be changed.



### NOTE:

Some of the BIOS options may vary throughout the life cycle of the product and are subject to change without prior notice.

---

### 5.1.1 Starting Setup

The UEFI BIOS is activated when the computer is turned on. The setup program can be activated in one of two ways.

1. Press the **DEL** key as soon as the system is turned on or
2. Press the **DEL** key when the “**Press DEL to enter SETUP**” message appears on the screen.

If the message disappears before the **DEL** key is pressed, restart the computer and try again.

### 5.1.2 Using Setup

Use the arrow keys to highlight items, press **ENTER** to select, use the PageUp and PageDown keys to change entries, press **F1** for help and press **Esc** to quit. Navigation keys are shown in the following table.

Key	Function
Up arrow	Move to the item above
Down arrow	Move to the item below
Left arrow	Move to the item on the left hand side
Right arrow	Move to the item on the right hand side
+	Increase the numeric value or make changes
-	Decrease the numeric value or make changes
F1 key	General help, only for Status Page Setup Menu and Option Page Setup Menu
F2 key	Load previous values.
F3 key	Load optimized defaults
F4 key	Save changes and Exit BIOS
Esc key	Main Menu – Quit and do not save changes into CMOS Status Page Setup Menu and Option Page Setup Menu -- Exit current page and return to Main Menu

**Table 5-1: BIOS Navigation Keys**

### 5.1.3 Getting Help

When **F1** is pressed a small help window describing the appropriate keys to use and the possible selections for the highlighted item appears. To exit the Help Window press **Esc** or the **F1** key again.

## POC-W22A-H81 Medical Panel PC

### 5.1.4 BIOS Menu Bar

The **menu bar** on top of the BIOS screen has the following main items:

- Main – Changes the basic system configuration.
- Advanced – Changes the advanced system settings.
- Chipset – Changes the chipset settings.
- Boot – Changes the system boot configuration.
- Security – Sets User and Supervisor Passwords.
- Save & Exit – Selects exit options and loads default settings
- Server Mgmt – Changes BMC network configuration.

The following sections completely describe the configuration options found in the menu items at the top of the BIOS screen and listed above.

## 5.2 Main

The **Main** BIOS menu (**BIOS Menu 1**) appears when the **BIOS Setup** program is entered. The **Main** menu gives an overview of the basic system information.

Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.						
Main	Advanced	Chipset	Boot	Security	Save & Exit	Server Mgmt
BIOS Information				Set the Date. Use Tab to switch between Data elements.		
BIOS Vendor	American Megatrends					
Core Version	4.6.5.4					
Compliance	UEFI 2.3.1; PI 1.2					
Project Version	E438AR32.ROM					
Build Date and Time	08/07/2015 11:25:07					
iWDD Vender		iEi				
iWDD Version		E438ER11.bin				
IPMI Module		N/A				
Processor Information				-----		
Name	Haswell		←→: Select Screen			
Brand String	Intel(R)Core(TM) i5-457		↑ ↓: Select Item			
Frequency	3200 MHz		EnterSelect			
Processor ID	306c3		+/-: Change Opt.			
Stepping	C0		F1: General Help			
Number of Processors	2Core(S)/4Thread(s)		F2: Previous Values			
Microcode Revision	1c		F3: Optimized Defaults			
GT Info	GT2 (700 MHz)		F4: Save & Exit			
IGF VBIOS Version		2179		ESC: Exit		
Memory RC Version		1.9.0.0				
Total Memory		4096 MB (DDR3)				
Memory Frequency		1600 Mhz				
PCH Information						
Name	LynxPoint					
PCH SKU	H81					
Stepping	05/C2					
LAN PHY Revision	A3					
ME FW Version		9.1.10.1005				
ME Firmware SKU		1.5MB				
SPI Clock Frequency						
D0FR Support	Unsupported					
Read Status Clock Frequency	50 MHz					
Write Status Clock Frequency	50 MHz					
Fast Read Status Clock Frequency	50 MHz					
System Date		[Tue 09/22/2015]				
System Time		[16:49:37]				
Access Level		Administrator				
Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.						

BIOS Menu 1: Main

## POC-W22A-H81 Medical Panel PC

### → System Date [xx/xx/xx]

Use the **System Date** option to set the system date. Manually enter the day, month and year.

### → System Time [xx:xx:xx]

Use the **System Time** option to set the system time. Manually enter the hours, minutes and seconds.

## 5.3 Advanced

Use the **Advanced** menu (**BIOS Menu 2**) to configure the CPU and peripheral devices through the following sub-menus:



### WARNING:

Setting the wrong values in the sections below may cause the system to malfunction. Make sure that the settings made are compatible with the hardware.

```

Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.
Main  Advanced  Chipset  Boot  Security  Save & Exit  Server Mgmt
-----
> ACPI Settings                System ACPI Parameters.
> RTC Wake Settings
> Trusted Computing
> CPU Configuration
> SATA Configuration
> USB Configuration
> F81866 Super IO Configuration
> iWDD H/M Monitor
> Serial Port Console Redirection
> iEi Feature

-----
<=>: Select Screen
↑ ↓: Select Item
Enter>Select
+/-: Change Opt.
F1:  General Help
F2:  Previous Values
F3:  Optimized Defaults
F4:  Save & Exit
ESC: Exit


Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.
    
```

### BIOS Menu 2: Advanced



### 5.3.1 ACPI Settings

The **ACPI Settings** menu (**BIOS Menu 3**) configures the Advanced Configuration and Power Interface (ACPI) options.



#### BIOS Menu 3: ACPI Settings

##### → ACPI Sleep State [S1 only (CPU Stop Clock)]


Use the **ACPI Sleep State** option to specify the sleep state the system enters when it is not being used.

- **S1 only (CPU Stop Clock) DEFAULT** The system enters S1 (POS) sleep state. The system appears off. The CPU is stopped; RAM is refreshed; the system is running in a low power mode.
- **S3 only (Suspend to RAM)** The caches are flushed and the CPU is powered off. Power to the RAM is maintained. The computer returns slower to a working state, but more power is saved.

## POC-W22A-H81 Medical Panel PC

### 5.3.2 RTC Wake Settings

The **RTC Wake Settings** menu (**BIOS Menu 4**) configures RTC wake event.



#### BIOS Menu 4: RTC Wake Settings

##### → Wake System with Fixed Time [Disabled]

Use the **Wake System with Fixed Time** option to specify the time the system should be roused from a suspended state.

- **Disabled**    **DEFAULT**    The real time clock (RTC) cannot generate a wake event

➔ **Enabled**


If selected, the following appears with values that can be selected:

- \*Wake up every day
- \*Wake up date
- \*Wake up hour
- \*Wake up minute
- \*Wake up second

After setting the alarm, the computer turns itself on from a suspend state when the alarm goes off.

### 5.3.3 Trusted Computing

Use the **Trusted Computing** menu (**BIOS Menu 5**) to configure settings related to the Trusted Computing Group (TCG) Trusted Platform Module (TPM).



**BIOS Menu 5: Trusted Computing**

## POC-W22A-H81 Medical Panel PC


### → Security Device Support [Enable]

Use the **Security Device Support** option to configure support for the security devices.

- **Disable**                      Security device support is disabled.
- **Enable**    **DEFAULT**      Security device support is enabled.

### 5.3.4 CPU Configuration

Use the **CPU Configuration (BIOS Menu 6)** to view detailed CPU specifications and configure the CPU.



Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.

Advanced

CPU Configuration		Enable for Windows XP and Linux (OS optimized for Hyper-Threading Technology) and Disabled for other OS (OS not optimized for Hyper-Threading Technology). When Disabled only one thread per enabled core is enabled.
Intel(R) Core(TM) i5-4570TE CPU @ 2.70GHz		
CPU Signature	306c3	
Microcode Patch	1c	
Max CPU Speed	2700 MHz	
Min CPU Speed	800 MHz	
CPU Speed	3200 MHz	
Processor Cores	2	
Intel HT Technology	Supported	
Intel VT-x Technology	Supported	
Intel SMX Technology	Supported	
64-bit	Supported	
EIST Technology	Supported	
-----		
L1 Data Cache	32 kB x 2	←→: Select Screen
L1 Code Cache	32 kB x 2	↑ ↓: Select Item
L2 Cache	256 kB x 2	Enter>Select
L3 Cache	4096 kB	+/-: Change Opt.
Hyper-threading	[Enabled]	F1: General Help
Active Processor cores	[All]	F2: Previous Values
Intel Virtualization Technology	[Disabled]	F3: Optimized Defaults
EIST	[Enabled]	F4: Save & Exit
		ESC: Exit

Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.

#### BIOS Menu 6: CPU Configuration

The CPU Configuration menu lists the following CPU details:

- Processor Type: Lists the brand name of the CPU being used
- CPU Signature: Lists the CPU signature value.
- Microcode Patch: Lists the microcode patch being used.

- Max CPU Speed: Lists the maximum CPU processing speed.
- Min CPU Speed: Lists the minimum CPU processing speed.
- CPU Speed: Lists the CPU processing speed.
- Processor Cores: Lists the number of the processor core
- Intel HT Technology: Indicates if Intel HT Technology is supported by the CPU.
- Intel VT-x Technology: Indicates if Intel VT-x Technology is supported by the CPU.
- Intel SMX Technology: Indicates if Intel SMX Technology is supported by the CPU.
- 64-bit: Indicates if 64-bit OS is supported by the CPU.
- EIST Technology: Indicates if EIST Technology is supported by the CPU.
- L1 Data Cache: Lists the amount of data storage space on the L1 cache.
- L1 Code Cache: Lists the amount of code storage space on the L1 cache.
- L2 Cache: Lists the amount of storage space on the L2 cache.
- L3 Cache: Lists the amount of storage space on the L3 cache.

➔ **Hyper-threading [Enabled]**

Use the **Hyper-threading** option to enable or disable the Intel® Hyper-Threading Technology.

- ➔ **Disabled** Disable Intel® Hyper-Threading Technology
- ➔ **Enabled** **DEFAULT** Enable Intel® Hyper-Threading Technology

➔ **Active Processor Cores [All]**

Use the **Active Processor Cores** BIOS option to enable numbers of cores in the processor package.

- ➔ **All** **DEFAULT** Enable all cores in the processor package.
- ➔ **1** Enable one core in the processor package.



## POC-W22A-H81 Medical Panel PC

### → Intel Virtualization Technology [Disabled]

Use the **Intel Virtualization Technology** option to enable or disable virtualization on the system. When combined with third party software, Intel® Virtualization technology allows several OSs to run on the same system at the same time.

- **Disabled**                      **DEFAULT**      Disables Intel Virtualization Technology.
- **Enabled**    Enables Intel Virtualization Technology.


### → EIST [Enabled]

Use the **EIST** option to enable or disable Enhanced Intel SpeedStep® Technology (EIST).

- **Disabled**    Disables Enhanced Intel SpeedStep® Technology.
- **Enabled**                      **DEFAULT**      Enables Enhanced Intel SpeedStep® Technology.

## 5.3.5 SATA Configuration

Use the **SATA Configuration** menu (**BIOS Menu 7**) to change and/or set the configuration of the SATA devices installed in the system.



**BIOS Menu 7: SATA Configuration**

➔ **SATA Controller(s) [Enabled]**

Use the **SATA Controller(s)** option to configure the SATA controller(s).

- ➔ **Enabled**      **DEFAULT**      Enable the on-board SATA controller(s).
- ➔ **Disabled**                      Disable the on-board SATA controller(s).


➔ **SATA Mode Selection [IDE]**

Use the **SATA Mode Selection** option to determine how SATA devices operate.

- ➔ **IDE**              **DEFAULT**      Configures SATA devices as normal IDE device.
- ➔ **AHCI**                      Configures SATA devices as AHCI device.

### 5.3.6 USB Configuration

Use the **USB Configuration** menu (**BIOS Menu 8**) to read USB configuration information and configure the USB settings.



**BIOS Menu 8: USB Configuration**

## POC-W22A-H81 Medical Panel PC

### → USB Devices

The **USB Devices Enabled** field lists the USB devices that are enabled on the system


### → Legacy USB Support [Enabled]

Use the **Legacy USB Support** BIOS option to enable USB mouse and USB keyboard support. Normally if this option is not enabled, any attached USB mouse or USB keyboard does not become available until a USB compatible operating system is fully booted with all USB drivers loaded. When this option is enabled, any attached USB mouse or USB keyboard can control the system even when there is no USB driver loaded onto the system.

- **Enabled**      **DEFAULT**      Legacy USB support enabled
- **Disabled**                      Legacy USB support disabled
- **Auto**                              Legacy USB support disabled if no USB devices are connected

### 5.3.7 F81866 Super IO Configuration


Use the **F81866 Super IO Configuration** menu (**BIOS Menu 9**) to set or change the configurations for the serial ports.



**BIOS Menu 9: F81866 Super IO Configuration**

### 5.3.7.1 Serial Port n Configuration

Use the **Serial Port n Configuration** menu (**BIOS Menu 10**) to configure the serial port n.



**BIOS Menu 10: Serial Port n Configuration Menu**

#### 5.3.7.1.1 Serial Port 1 Configuration

➔ **Serial Port [Enabled]**

Use the **Serial Port** option to enable or disable the serial port.

- ➔ **Disabled**                      Disable the serial port
- ➔ **Enabled      DEFAULT**      Enable the serial port

➔ **Change Settings [Auto]**

Use the **Change Settings** option to change the serial port IO port address and interrupt address.

- ➔ **Auto                      DEFAULT**      The serial port IO port address and interrupt address are automatically detected.
- ➔ **IO=3F8h; IRQ=4**                      Serial Port I/O port address is 3F8h and the interrupt address is IRQ4

## POC-W22A-H81 Medical Panel PC

- **IO=3F8h;**  
**IRQ=3, 4**                      Serial Port I/O port address is 3F8h and the interrupt address is IRQ3, 4
- **IO=2C0h;**  
**IRQ=3, 4**                      Serial Port I/O port address is 2C0h and the interrupt address is IRQ3, 4
- **IO=2C8h;**  
**IRQ=3, 4**                      Serial Port I/O port address is 2C8h and the interrupt address is IRQ3, 4

### → **Device Mode [RS232]**

Use the **Device Mode** option to set the Serial Port 1 signaling mode.

- **RS232**            **DEFAULT**            Configure Serial Port 1 as RS-232
- **RS422/485**                      Configure Serial Port 1 as RS-422/485

### 5.3.7.1.2 Serial Port 2 Configuration

#### → **Serial Port [Enabled]**

Use the **Serial Port** option to enable or disable the serial port.

- **Disabled**                      Disable the serial port
- **Enabled**            **DEFAULT**            Enable the serial port

#### → **Change Settings [Auto]**

Use the **Change Settings** option to change the serial port IO port address and interrupt address.

- **Auto**                      **DEFAULT**            The serial port IO port address and interrupt address are automatically detected.
- **IO=2F8h; IRQ=3**                      Serial Port I/O port address is 2F8h and the interrupt address is IRQ3
- **IO=3F8h;**  
**IRQ=3, 4**                      Serial Port I/O port address is 3F8h and the interrupt address is IRQ3, 4



- ➔ **IO=2F8h;**  
**IRQ=3, 4**  
Serial Port I/O port address is 2F8h and the interrupt address is IRQ3, 4
- ➔ **IO=2C0h;**  
**IRQ=3, 4**  
Serial Port I/O port address is 2C0h and the interrupt address is IRQ3, 4
- ➔ **IO=2C8h;**  
**IRQ=3, 4**  
Serial Port I/O port address is 2C8h and the interrupt address is IRQ3, 4

### 5.3.8 iWDD H/W Monitor

The **iWDD H/W Monitor** menu (**BIOS Menu 11**) shows the operating temperatures and voltages.

```
Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.
Advanced
PC Health Status
Smart Fan Mode Select
>Smart Fan Mode Configuration
CPU temperature           : +62 °C
System temperature       : +48 °C
CPU_FAN1 Speed           : 1108 RPM
SYS_FAN1 Speed           : N/A
CPU_CORE                 : +1.705 V
+5V                      : +4.953 V
+12V                    : +11.866 V
DDR                      : +1.496 V
+5VSB                   : +6.120 V
+3.3V                   : +3.230 V
+3.3VSB                 : +3.207 V
-----
<->: Select Screen
↑ ↓: Select Item
Enter>Select
+/-: Change Opt.
F1:  General Help
F2:  Previous Values
F3:  Optimized Defaults
F4:  Save & Exit
ESC: Exit
Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.
```

#### BIOS Menu 11: iWDD H/W Monitor

#### ➔ Hardware Health Status

The following system parameters and values are shown. The system parameters that are monitored are:

- Temperature:
  - CPU Temperature
  - System Temperature

## POC-W22A-H81 Medical Panel PC

- Fan Speed
  - CPU fan
  - System fan
- Voltages:
  - CPU\_CORE
  - +5V
  - +12V
  - +DDR
  - +5VSB
  - +3.3V
  - +3.3VSB

### 5.3.8.1 Smart Fan Mode Configuration

Use the **Smart Fan Mode Configuration** submenu (**BIOS Menu 12**) to configure fan speed settings.

```

Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.
  Advanced
Smart Fan Mode Configuration
CPU_FAN1 Smart Fan Control      [Auto Mode]
Fan start temperature           50
Fan off temperature             40
Fan start PWM                   30
Fan slope PWM                   1
SYS_FAN1 Smart Fan Control      [Auto Mode]
Fan start temperature           50
Fan off temperature             40
Fan start PWM                   30
Fan slope PWM                   1
Smart Fan Mode Select
-----
<->: Select Screen
↑ ↓: Select Item
Enter: Select
+/-: Change Opt.
F1: General Help
F2: Previous Values
F3: Optimized Defaults
F4: Save & Exit
ESC: Exit
Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.
  
```

#### BIOS Menu 12: Smart Fan Mode Configuration

➔ **CPU\_FAN1 Smart Fan Control/SYS\_FAN1 Smart Fan Control [Auto Mode]**

Use the **CPU\_FAN1 Smart Fan Control/SYS\_FAN1 Smart Fan Control** option to configure the CPU/System Smart Fan.

- **Manual Mode**                      The fan spins at the speed set in Manual Mode settings.
- **Auto Mode**                      **DEFAULT**      The fan adjusts its speed using Auto Mode settings.

→ **Fan start/off temperature**

Use the + or – key to change the **Fan start/off temperature** value. Enter a decimal number between 1 and 100.

→ **Fan start PWM**

Use the + or – key to change the **Fan start PWM** value. Enter a decimal number between 1 and 100.

→ **Fan slope PWM**

Use the + or – key to change the **Fan slope PWM** value. Enter a decimal number between 1 and 8.

## POC-W22A-H81 Medical Panel PC

### 5.3.9 Serial Port Console Redirection

The **Serial Port Console Redirection** menu (**BIOS Menu 13**) allows the console redirection options to be configured. Console redirection allows users to maintain a system remotely by re-directing keyboard input and text output through the serial port.

```

Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.
  Advanced
COM1
Console Redirection          [Disabled]      Console Redirection
> Console Redirection Settings      Enable or Disable.

COM2
Console Redirection          [Disabled]
> Console Redirection Settings
-----
<->: Select Screen
↑ ↓: Select Item
Enter>Select
+/-: Change Opt.
F1:  General Help
F2:  Previous Values
F3:  Optimized Defaults
F4:  Save & Exit
ESC: Exit

COM3 (BMC) (Disabled)
Console Redirection          Port Is Disabled

Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.
  
```

#### BIOS Menu 13: Serial Port Console Redirection


##### → Console Redirection [Disabled]

Use **Console Redirection** option to enable or disable the console redirection function.

- **Disabled**     **DEFAULT**     Disabled the console redirection function
- **Enabled**                     Enabled the console redirection function

### 5.3.10 IEI Feature

Use the **IEI Feature** menu (**BIOS Menu 14**) to configure One Key Recovery function.



#### BIOS Menu 14: IEI Feature

##### ➔ Auto Recovery Function [Disabled]

Use the **Auto Recovery Function** BIOS option to enable or disable the auto recovery function of the IEI One Key Recovery.

- ➔ **Disabled**     **DEFAULT**     Auto recovery function disabled
- ➔ **Enabled**     Auto recovery function enabled



## POC-W22A-H81 Medical Panel PC

### 5.4 Chipset

Use the **Chipset** menu (**BIOS Menu 15**) to configure the system chipset.

```

Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.
Main   Advanced  Chipset   Boot   Security  Save & Exit  Server Mgmt
-----
> PCH-IO Configuration
> System Agent (SA) Configuration

PCH Parameters
-----
<->: Select Screen
↑ ↓: Select Item
Enter>Select
+ - Change Opt.
F1  General Help
F2  Previous Values
F3  Optimized Defaults
F4  Save & Exit
ESC Exit

Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.
  
```

**BIOS Menu 15: Chipset**

#### 5.4.1 PCH-IO Configuration

Use the **PCH-IO Configuration** menu (**BIOS Menu 16**) to configure the PCH-IO chipset.

```

Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.
-----
Chipset
-----
Auto Power Button Status      [Disable (ATX)]
Restore AC Power Loss         [Last State]

> PCI Express Configuration
> PCH Azalia Configuration

Select AC power state when
power is re-applied after
a power failure.
-----
<->: Select Screen
↑ ↓: Select Item
Enter>Select
+ - Change Opt.
F1  General Help
F2  Previous Values
F3  Optimized Defaults
F4  Save & Exit
ESC Exit

Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.
  
```

**BIOS Menu 16: PCH-IO Configuration**


➔ **Restore AC Power Loss [Last State]**

Use the **Restore AC Power** BIOS option to specify what state the system returns to if there is a sudden loss of power to the system.

- ➔ **Power Off**                      The system remains turned off
- ➔ **Power On**                        The system turns on
- ➔ **Last State**    **DEFAULT**    The system returns to its previous state. If it was on, it turns itself on. If it was off, it remains off.

### 5.4.1.1 PCI Express Configuration

Use the **PCI Express Configuration** submenu (**BIOS Menu 17**) to configure the PCI Express slots.



**BIOS Menu 17: PCI Express Configuration**

The **M\_PCIE\_1 Slot** and **M\_PCIE\_2 Slot** submenus both contain the following options:

➔ **PCIe Speed [Auto]**

Use the **PCIe Speed** option to configure the PCIe interface speed.

- Auto                      **DEFAULT**
- Gen 2

## POC-W22A-H81 Medical Panel PC

- Gen 1


### → Detect Non-Compliance Device [Disabled]

Use the **Detect Non-Compliance Device** option to enable or disable detecting if a non-compliance PCI Express device is connected to the PCI Express slot.

- **Disabled**                      **DEFAULT**      Disables to detect if a non-compliance PCI Express device is connected to the PCI Express slot.
- **Enabled**                                      Enables to detect if a non-compliance PCI Express device is connected to the PCI Express slot.

### 5.4.1.2 PCH Azalia Configuration

Use the **PCH Azalia Configuration** submenu (**BIOS Menu 18**) to configure the High Definition Audio codec.



**BIOS Menu 18: PCH Azalia Configuration**


➔ **Azalia(HD Audio) [Enabled]**

Use the **Azalia(HD Audio)** BIOS option to enable or disable the High Definition Audio controller.

- ➔ **Disabled**                      The High Definition Audio controller is disabled.
- ➔ **Enabled      DEFAULT**      The High Definition Audio controller is enabled.

### 5.4.2 System Agent (SA) Configuration

Use the **System Agent (SA) Configuration** menu (**BIOS Menu 19**) to configure the System Agent (SA) parameters.



#### BIOS Menu 19: System Agent (SA) Configuration

➔ **VT-d [Disabled]**

Use the **VT-d** option to enable or disable VT-d support.

- ➔ **Disabled      DEFAULT**      Disable VT-d support.
- ➔ **Enabled**                      Enable VT-d support.

## POC-W22A-H81 Medical Panel PC

### 5.4.2.1 Graphics Configuration

Use the **Graphics Configuration** menu (**BIOS Menu 20**) to configure the graphics settings.

```

Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.
Chipset
Graphics Configuration
Primary Display                [Auto]
DVMT Pre-Allocated            [256M]
DVMT Total Gfx Mem            [MAX]
> LCD Control

Select which of
Auto/IGFX/PCIE Graphics
device should be Primary
Display Or select SG for
Switchable Gfx.

-----
<=>: Select Screen
↑↓: Select Item
Enter>Select
+ - Change Opt.
F1  General Help
F2  Previous Values
F3  Optimized Defaults
F4  Save & Exit
ESC Exit

Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.
    
```

#### BIOS Menu 20: Graphics Configuration

##### → Primary Display [Auto]

Use the **Primary Display** option to select the graphics controller used as the primary boot device. Configuration options are listed below:

- Auto                   **DEFAULT**
- IGFX
- PCIE

##### → DVMT Pre-Allocated [256M]

Use the **DVMT Pre-Allocated** option to set the amount of system memory allocated to the integrated graphics processor when the system boots. The system memory allocated can then only be used as graphics memory, and is no longer available to applications or the operating system. Configuration options are listed below:

- 32M



- 64M
- 128M
- 256M           **DEFAULT**
- 512M

→ **DVMT Total Gfx Mem [MAX]**

Use the **DVMT Total Gfx Mem** option to select DVMT5.0 total graphic memory size used by the internal graphic device. The following options are available:

- 128M
- 256M
- MAX           **DEFAULT**

→ **Primary IGFX Boot Display [VBIOS Default]**

Use the **Primary IGFX Boot Display** option to select the display device used by the system when it boots. Configuration options are listed below.

- VBIOS Default   **DEFAULT**
- CRT
- LVDS
- HDMI

## POC-W22A-H81 Medical Panel PC

### 5.4.2.2 Memory Configuration

Use the **Memory Configuration** submenu (**BIOS Menu 21**) to display the memory information.

```

Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.
  Chipset
-----
Memory Information
Total Memory          4096 MB (DDR3)
CHA_DIMM 1           2048 MB (DDR3)
CHB_DIMM 1           2048 MB (DDR3)
-----
<->: Select Screen
↑ ↓: Select Item
EnterSelect
+ - Change Opt.
F1  General Help
F2  Previous Values
F3  Optimized Defaults
F4  Save & Exit
ESC Exit

Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.
  
```

**BIOS Menu 21: Memory Configuration**

## 5.5 Boot

Use the **Boot** menu (**BIOS Menu 22**) to configure system boot options.

```

Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.
  Main  Advanced  Chipset  Boot  Security  Save & Exit  Server Mgmt
-----
Boot Configuration
Bootup NumLock State      [On]
Quiet Boot                 [Enabled]
Option ROM Messages       [Force BIOS]
Launch PXE OpROM          [Disabled]
UEFI Boot                  [Disabled]
-----
Boot Option Priorities
-----
Select the keyboard
NumLock state
-----
<->: Select Screen
↑ ↓: Select Item
EnterSelect
+/-: Change Opt.
F1:  General Help
F2:  Previous Values
F3:  Optimized Defaults
F4:  Save & Exit
ESC: Exit

Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.
  
```

**BIOS Menu 22: Boot**

➔ **Bootup NumLock State [On]**

Use the **Bootup NumLock State** BIOS option to specify if the number lock setting must be modified during boot up.

➔ **On**                      **DEFAULT**                      Allows the Number Lock on the keyboard to be enabled automatically when the computer system boots up. This allows the immediate use of the 10-key numeric keypad located on the right side of the keyboard. To confirm this, the Number Lock LED light on the keyboard is lit.

➔ **Off**    Does not enable the keyboard Number Lock automatically. To use the 10-keys on the keyboard, press the Number Lock key located on the upper left-hand corner of the 10-key pad. The Number Lock LED on the keyboard lights up when the Number Lock is engaged.

➔ **Quiet Boot [Enabled]**

Use the **Quiet Boot** BIOS option to select the screen display when the system boots.

➔ **Disabled**                                      Normal POST messages displayed

➔ **Enabled**                      **DEFAULT**                      OEM Logo displayed instead of POST messages

➔ **Option ROM Messages [Force BIOS]**

Use the **Option ROM Messages** option to set the Option ROM display mode.

➔ **Force BIOS**                      **DEFAULT**                      Sets display mode to force BIOS.

➔ **Keep Current**                                      Sets display mode to current.

## POC-W22A-H81 Medical Panel PC

### → Launch PXE OpROM [Disabled]

Use the **Launch PXE OpROM** option to enable or disable boot option for legacy network devices.

- **Disabled**      **DEFAULT**      Ignore all PXE Option ROMs
- **Enabled**                      Load PXE Option ROMs

### → UEFI Boot [Disabled]

Use the **UEFI Boot** BIOS option to enable or disable UEFI boot.

- **Disabled**      **DEFAULT**      Disable UEFI boot.
- **Enabled**                      Enable UEFI boot if the 1<sup>st</sup> boot device is a GPT HDD.

## 5.6 Security

Use the **Security** menu (**BIOS Menu 23**) to set system and user passwords.

```

Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.
Main      Advanced  Chipset   Boot     Security  Save & Exit  Server Mgmt
-----
Password Description
If ONLY the Administrator's password is set,
then this only limits access to Setup and is
only asked for when entering Setup.
If ONLY the User's password is set, then this
is a power on password and must be entered to
boot or enter Setup. In Setup the User will
have Administrator rights.
The password must be
In the following range:
Maximum length                3
Minimum length                20

Administrator Password
User Password

Set Administrator
Password

-----
<->: Select Screen
↑ ↓: Select Item
EnterSelect
+/-: Change Opt.
F1: General Help
F2: Previous Values
F3: Optimized Defaults
F4: Save & Exit
ESC: Exit

Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.
    
```

**BIOS Menu 23: Security**

➔ **Administrator Password**


Use the **Administrator Password** field to set or change an administrator password.

➔ **User Password**

Use the **User Password** field to set or change a user password.

## 5.7 Save & Exit

Use the **Save & Exit** menu (**BIOS Menu 24**) to load default BIOS values, optimal failsafe values and to save configuration changes.



**BIOS Menu 24: Save & Exit**

➔ **Save Changes and Reset**

Use the **Save Changes and Reset** option to save the changes made to the BIOS options and reset the system.

➔ **Discard Changes and Reset**

Use the **Discard Changes and Reset** option to exit the system without saving the changes made to the BIOS configuration setup program.



## POC-W22A-H81 Medical Panel PC

### → Restore Defaults

Use the **Restore Defaults** option to load the optimal default values for each of the parameters on the Setup menus. **F3 key can be used for this operation.**

### → Save as User Defaults

Use the **Save as User Defaults** option to save the changes done so far as user defaults.

### → Restore User Defaults

Use the **Restore User Defaults** option to restore the user defaults to all the setup options.

## 5.8 Server Management

Use the **Server Mgmt** menu (**BIOS Menu 25**) to configure BMC network parameters.

```

Aptio Setup Utility - Copyright (C) 2015 American Megatrends, Inc.
Main   Advanced  Chipset  Boot   Security  Save & Exit  Server Mgmt
-----
BMC Self Test Status          FAILED
> System Event Log
> BMC network configuration

Press <Enter> to change
the SEL event log
configuration.


-----
<=>: Select Screen
↑ ↓: Select Item
EnterSelect
+ - Change Opt.
F1  General Help
F2  Previous Values
F3  Optimized Defaults
F4  Save & Exit
ESC Exit

Version 2.17.1246. Copyright (C) 2015 American Megatrends, Inc.
    
```

**BIOS Menu 25: Server Management**

### 5.8.1 System Event Log

Use the **System Event Log** menu (**BIOS Menu 26**) to configure the system event log of the BMC.



#### BIOS Menu 26: System Event Log

##### → SEL Components [Enabled]

Use the **SEL Components** option to enable or disable all features of System Event Log (SEL) when the system boots.

- **Disabled** System Event Log is disabled.
- **Enabled** **DEFAULT** System Event Log is enabled.

##### → Erase SEL [No]

Use the **Erase SEL** option for erasing SEL.

- **No** **DEFAULT** Do not erase system event log.
- **Yes, On next reset** Erase system event log on next reset.

## POC-W22A-H81 Medical Panel PC

→ **Yes, On every reset** Erase system event log on every reset.

→ **When SEL is Full [Do Nothing]**

Use the **When SEL is Full** option to select an reaction to a full SEL.

→ **Do Nothing**      **DEFAULT**      Do not do anything when SEL is full.

→ **Erase Immediately** Erase SEL immediately when SEL is full.

→ **Log EFI Status Codes [Both]**

Use the **Log EFI Status Codes** option to configure how to log Extensible Firmware Interface (EFI) status codes.

→ **Disabled** Disable the logging of EFI status codes.


→ **Both**      **DEFAULT**      Log both error codes and progress codes.

→ **Error code** Log error codes only.

→ **Progress code** Log progress codes only.

### 5.8.2 BMC Network Configuration

Use the **BMC Network Configuration** menu (**BIOS Menu 27**) to configure the BMC network parameters.



#### BIOS Menu 27: BMC Network Configuration

##### → Configuration Address source [Unspecified]

Use the **Configuration Address source** option to configure LAN channel parameters statically or dynamically (by BIOS or BMC). Choosing the **Unspecified** option will not modify any BMC network parameters during BIOS phase. The following options are available:

- Unspecified     **DEFAULT**
- Static
- Dynamic-Obtained by BMC
- Dynamic-Loaded by BIOS
- Dynamic-BMC running other Protocol

Chapter

**6**

# System Maintenance

---



## 6.1 System Maintenance Introduction

If the components of the POC-W22A-H81 fail they must be replaced. Please contact the system reseller or vendor to purchase the replacement parts.

## 6.2 Anti-static Precautions

---



### **WARNING:**

Failure to take ESD precautions during the maintenance of the POC-W22A-H81 may result in permanent damage to the POC-W22A-H81 and severe injury to the user.

---

Electrostatic discharge (ESD) can cause serious damage to electronic components, including the POC-W22A-H81. Dry climates are especially susceptible to ESD. It is therefore critical that whenever the POC-W22A-H81 is accessed internally, or any other electrical component is handled, the following anti-static precautions are strictly adhered to.

- ***Wear an anti-static wristband:*** - Wearing a simple anti-static wristband can help to prevent ESD from damaging the board.
- ***Self-grounding:*** - Before handling the board touch any grounded conducting material. During the time the board is handled, frequently touch any conducting materials that are connected to the ground.
- ***Use an anti-static pad:*** - When configuring the POC-W22A-H81, place it on an anti-static pad. This reduces the possibility of ESD damaging the POC-W22A-H81.
- ***Only handle the edges of the PCB:*** - When handling the PCB, hold the PCB by the edges.

## POC-W22A-H81 Medical Panel PC

### 6.3 Turn off the Power



#### **WARNING:**

Failing to turn off the system before opening it can cause permanent damage to the system and serious or fatal injury to the user.

Before any maintenance procedures are carried out on the system, make sure the system is turned off.

### 6.4 System Fan Replacement

The POC-W22A-H81 has one system fan which can be accessed by removing the back cover. To replace the system fan, follow the instructions below.

- Step 1:** Follow all anti-static procedures. See **Section 6.2**.
- Step 2:** Turn off the power. See **Section 6.3**.
- Step 3:** Remove the back cover. See **Section 3.4**.
- Step 4:** Remove the three system fan retention screws (M3\*8) and disconnect the fan cable. Lift the fan from the system.




Figure 6-1: System Fan Retention Screws and Connector

**Step 5:** Install a new system fan and secure it with the previously removed retention screws.

**Step 6:** Re-install the back cover.

## 6.5 SO-DIMM Replacement


The POC-W22A-H81 has two SO-DIMM slots. Each SO-DIMM slot is pre-installed with one 2 GB DDR3 SO-DIMM. To replace the SO-DIMM, follow the instructions below.

**Step 1:** Follow all anti-static procedures. See **Section 6.2**.

**Step 2:** Turn off the power. See **Section 6.3**.

**Step 3:** Remove the back cover. See **Section 3.4**.

**Step 4:** Locate the SO-DIMM slot cover (**Figure 6-2**).



**Figure 6-2: SO-DIMM Slot Cover Location**

**Step 5:** To access the SO-DIMMs, simply open the SO-DIMM slot cover (See **Figure 6-3**).

## POC-W22A-H81 Medical Panel PC




Figure 6-3: SO-DIMM Locations

- Step 6:** Remove the SO-DIMM by releasing the arms on the SO-DIMM socket.
- Step 7:** Align the new SO-DIMM with the socket. The SO-DIMM must be oriented in such a way that the notch in the middle of the SO-DIMM must be aligned with the plastic bridge in the socket (**Figure 6-4**).
- Step 8:** Insert the SO-DIMM. Push the SO-DIMM chip into the socket at an angle (**Figure 6-4**).




Figure 6-4: SO-DIMM Installation

- Step 9:** Open the SO-DIMM socket arms. Gently pull the arms of the SO-DIMM socket out and push the rear of the SO-DIMM down (See **Figure 6-4**).



**Step 10:** Secure the SO-DIMM. Release the arms on the SO-DIMM socket. They clip into place and secure the SO-DIMM in the socket.

**Step 11:** Close the SO-DIMM slot cover.

**Step 12:** Re-install the back cover.

## 6.6 Reinstalling the Back Cover

---



### **WARNING:**

Failing to reinstall the back cover may result in permanent damage to the system. Please make sure all coverings are properly installed.

---

When maintenance procedures are complete, please make sure the plastic back cover is replaced.

## 6.7 Troubleshooting

If the following situations happen, contact your distributor, sales representatives or IEI customer service center for technical support.

- The HDD is installed correctly, but the POC-W22A-H81 is unable to boot with AC power input after pressing the power button.
- Unable to shut down the POC-W22A-H81 normally
- The system fan makes a loud and annoying noise

Please have the following information prepared prior to reporting the abnormal situations:

- Product name and S/N
- OS, BIOS version and applications
- A complete description of the abnormal situation (with photos or video if available)



Appendix

**A**

# Regulatory Compliance

---

**DECLARATION OF CONFORMITY**

This equipment is in conformity with the following EU directives:

- EMC Directive (2004/108/EC, 2014/30/EU)
- Low-Voltage Directive (2006/95/EC, 2014/35/EU)
- RoHS II Directive (2011/65/EU, 2015/863/EU)
- Medical Device Directive 93/42/EEC: EN 60601-1

If the user modifies and/or install other devices in the equipment, the CE conformity declaration may no longer apply.

If this equipment has telecommunications functionality, it also complies with the requirements of the Radio Equipment Directive 2014/53/EU.

---

**English**

IEI Integration Corp declares that this equipment is in compliance with the essential requirements and other relevant provisions of Directive 2014/53/EU.

---

**Български [Bulgarian]**

IEI Integration Corp. декларира, че този оборудване е в съответствие със съществените изисквания и другите приложими правила на Директива 2014/53/EU.

---

**Česky [Czech]**

IEI Integration Corp tímto prohlašuje, že tento zařizení je ve shodě se základními požadavky a dalšími příslušnými ustanoveními směrnice 2014/53/EU.

---

**Dansk [Danish]**

IEI Integration Corp erklærer herved, at følgende udstyr overholder de væsentlige krav og øvrige relevante krav i direktiv 2014/53/EU.

---

**Deutsch [German]**

IEI Integration Corp, erklärt dieses Gerät entspricht den grundlegenden Anforderungen und den weiteren entsprechenden Vorgaben der Richtlinie 2014/53/EU.

---

**Eesti [Estonian]**

IEI Integration Corp deklareerib seadme seadme vastavust direktiivi 2014/53/EÜ põhinõuetele ja nimetatud direktiivist tulenevatele teistele asjakohastele sätetele.

---

## POC-W22A-H81 Medical Panel PC

---

### Español [Spanish]

IEI Integration Corp declara que el equipo cumple con los requisitos esenciales y cualesquiera otras disposiciones aplicables o exigibles de la Directiva 2014/53/EU.

---

### Ελληνική [Greek]

IEI Integration Corp ΔΗΛΩΝΕΙ ΟΤΙ ΕΞΟΠΛΙΣΜΟΣ ΣΥΜΜΟΡΦΩΝΕΤΑΙ ΠΡΟΣ ΤΙΣ ΟΥΣΙΩΔΕΙΣ ΑΠΑΙΤΗΣΕΙΣ ΚΑΙ ΤΙΣ ΛΟΙΠΕΣ ΣΧΕΤΙΚΕΣ ΔΙΑΤΑΞΕΙΣ ΤΗΣ ΟΔΗΓΙΑΣ 2014/53/EU.

---

### Français [French]

IEI Integration Corp déclare que l'appareil est conforme aux exigences essentielles et aux autres dispositions pertinentes de la directive 2014/53/EU.

---

### Italiano [Italian]

IEI Integration Corp dichiara che questo apparecchio è conforme ai requisiti essenziali ed alle altre disposizioni pertinenti stabilite dalla direttiva 2014/53/EU.

---

### Latviski [Latvian]

IEI Integration Corp deklarē, ka iekārta atbilst būtiskajām prasībām un citiem ar to saistītajiem noteikumiem Direktīvas 2014/53/EU.

---

### Lietuvių [Lithuanian]

IEI Integration Corp deklaruoja, kad šis įranga atitinka esminius reikalavimus ir kitas 2014/53/EU Direktyvos nuostatas.

---

### Nederlands [Dutch]

IEI Integration Corp dat het toestel toestel in overeenstemming is met de essentiële eisen en de andere relevante bepalingen van richtlijn 2014/53/EU.

---

### Malti [Maltese]

IEI Integration Corp jiddikjara li dan prodott jikkonforma mal-ħtiġijiet essenzjali u ma provvedimenti oħrajn relevanti li hemm fid-Dirrettiva 2014/53/EU.

---

### Magyar [Hungarian]

IEI Integration Corp nyilatkozom, hogy a berendezés megfelel a vonatkozó alapvető követelményeknek és az 2014/53/EU irányelv egyéb előírásainak.

---

### Polski [Polish]

IEI Integration Corp oświadcza, że wyrobu jest zgodny z zasadniczymi wymogami oraz pozostałymi stosownymi postanowieniami Dyrektywy 2014/53/EU.

---

### Português [Portuguese]

IEI Integration Corp declara que este equipamento está conforme com os requisitos essenciais e outras disposições da Directiva 2014/53/EU.

---

---

**Româna [Romanian]**

IEI Integration Corp declară că acest echipament este în conformitate cu cerințele esențiale și cu celelalte prevederi relevante ale Directivei 2014/53/EU.

---

**Slovensko [Slovenian]**

IEI Integration Corp izjavlja, da je ta opreme v skladu z bistvenimi zahtevami in ostalimi relevantnimi določili direktive 2014/53/EU.

---

**Slovensky [Slovak]**

IEI Integration Corp týmto vyhlasuje, že zariadenia spĺňa základné požiadavky a všetky príslušné ustanovenia Smernice 2014/53/EU.

---

**Suomi [Finnish]**

IEI Integration Corp vakuuttaa täten että laitteet on direktiivin 2014/53/EU oleellisten vaatimusten ja sitä koskevien direktiivin muiden ehtojen mukainen.

---

**Svenska [Swedish]**

IEI Integration Corp förklarar att denna utrustningstyp står i överensstämmelse med de väsentliga egenskapskrav och övriga relevanta bestämmelser som framgår av direktiv 2014/53/EU.

---

**FCC WARNING**

This equipment complies with part 18 of the FCC Rules.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

## POC-W22A-H81 Medical Panel PC

### UL CLASSIFIED



Medical general medical equipment with respect to electrical shock, fire and mechanical hazards only in accordance with ANSI/AAMI ES60601-1 (2005 and Amendment 1), CAN/CSA-C22.2 NO. 60601-1 (2014).

### ROHS STATEMENT



The label on the product indicates this product conforms to European (EU) Restriction of Hazardous Substances (RoHS) that set maximum concentration limits on hazardous materials used in electrical and electronic equipment.

### CHINA ROHS



The label on the product indicates the estimated “Environmentally Friendly Use Period” (EFUP). This is an estimate of the number of years that these substances would “not leak out or undergo abrupt change.” This product may contain replaceable sub-assemblies/components which have a shorter EFUP such as batteries and lamps. These components will be separately marked.



Appendix

**B**

# Safety Precautions

---

**WARNING:**

The precautions outlined in this chapter should be strictly followed. Failure to follow these precautions may result in permanent damage to the POC-W22A-H81.

## B.1 Safety Precautions

Please follow the safety precautions outlined in the sections that follow:

### B.1.1 General Safety Precautions

Please ensure the following safety precautions are adhered to at all times.

- ***To prevent the risk of electric shock, make sure power cord is unplugged from wall socket.*** To fully disengage the power to the unit, please disconnect the power cord from the ac outlet. Refer servicing to qualified service personnel. The AC outlet shall be readily available and accessible.
- ***Users must not allow SIP/SOPs and the patient to come into contact at the same time.***
- ***Grounding reliability*** can only be achieved when the equipment is connected to an equivalent receptacle marked “Hospital Only” or “Hospital Grade”.
- ***Follow the electrostatic precautions*** outlined below whenever the POC-W22A-H81 is opened.
- ***Make sure the power is turned off and the power cord is disconnected*** whenever the POC-W22A-H81 is being installed, moved or modified.
- ***Do not apply voltage levels that exceed the specified voltage range.*** Doing so may cause fire and/or an electrical shock. Use a power cord that matches the voltage of the power outlet, which has been approved and complies with the safety standard of your particular country.

- **Electric shocks can occur** if the POC-W22A-H81 chassis is opened when the POC-W22A-H81 is running. To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- **Do not drop or insert any objects** into the ventilation openings of the POC-W22A-H81.
- **If considerable amounts of dust, water, or fluids enter the POC-W22A-H81**, turn off the power supply immediately, unplug the power cord, and contact the POC-W22A-H81 vendor.
- **DO NOT:**
  - Drop the POC-W22A-H81 against a hard surface.
  - Strike or exert excessive force onto the LCD panel.
  - Touch any of the LCD panels with a sharp object
  - In a site where the ambient temperature exceeds the rated temperature

### B.1.2 Anti-static Precautions

---



#### **WARNING:**

Failure to take ESD precautions during the installation of the POC-W22A-H81 may result in permanent damage to the POC-W22A-H81 and severe injury to the user.

---

Electrostatic discharge (ESD) can cause serious damage to electronic components, including the POC-W22A-H81. Dry climates are especially susceptible to ESD. It is therefore critical that whenever the POC-W22A-H81 is opened and any of the electrical components are handled, the following anti-static precautions are strictly adhered to.

- **Wear an anti-static wristband:** Wearing a simple anti-static wristband can help to prevent ESD from damaging any electrical component.
- **Self-grounding:** Before handling any electrical component, touch any grounded conducting material. During the time the electrical component is handled, frequently touch any conducting materials that are connected to the ground.

## POC-W22A-H81 Medical Panel PC

- **Use an anti-static pad:** When configuring or working with an electrical component, place it on an anti-static pad. This reduces the possibility of ESD damage.
- **Only handle the edges of the electrical component:** When handling the electrical component, hold the electrical component by its edges.

### B.1.3 Product Disposal

---



#### CAUTION:

Risk of explosion if battery is replaced by an incorrect type. Only certified engineers should replace the on-board battery.

Dispose of used batteries according to instructions and local regulations.

---

- Outside the European Union - If you wish to dispose of used electrical and electronic products outside the European Union, please contact your local authority so as to comply with the correct disposal method.
- Within the European Union - The device that produces less waste and is easier to recycle is classified as electronic device in terms of the European Directive 2012/19/EU (WEEE), and must not be disposed of as domestic garbage.



EU-wide legislation, as implemented in each Member State, requires that waste electrical and electronic products carrying the mark (left) must be disposed of separately from normal household waste. This includes monitors and electrical accessories, such as signal cables or power cords.

When you need to dispose of your display products, please follow the guidance of your local authority, or ask the shop where you purchased the product. The mark on electrical and electronic products only applies to the current European Union Member States.

Please follow the national guidelines for electrical and electronic product disposal.

### B.1.4 Classification

- Power by Class I power supply (IEI, POC-W22A-H81)
- No Applied Part.
- No protection against the ingress of water: IPX0
- Mode of operation: Continuous Operation

The equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide: Not AP or APG Category.

## B.2 Maintenance and Cleaning Precautions

When maintaining or cleaning the POC-W22A-H81, please follow the guidelines below.



### **WARNING:**

If you dropped any material or liquid such as water onto the panel PC when cleaning, unplug the power cable immediately and contact your dealer or the nearest service center. Always make sure your hands are dry when unplugging the power cable.



### **CAUTION:**

- For safety reasons, turn-off the power switch and unplug the panel PC before cleaning.
- Do not scratch or rub the screen with a hard object.
- Never use any of the following solvents on the medical panel PC. Harsh chemicals may cause damage to the cabinet and the touch sensor.

**Thinner Spray-type cleaner, Benzene, Wax, Abrasive cleaner, Acid or Alkaline solvent.**

---

### B.2.1 Maintenance and Cleaning

Prior to cleaning any part or component of the POC-W22A-H81, please read the details below.



## POC-W22A-H81 Medical Panel PC

- To clean the POC-W22A-H81,
  - remove dirt with a lightly moistened cloth. Then wipe the external chassis with a soft dry cloth.
  - use 75% ethanol alcohol to clean the external chassis.
- Cleaning frequency: follow the cleaning method guidelines of the hospital.
- Except for the LCD panel, never spray or squirt liquids directly onto any other components.
- The interior of the POC-W22A-H81 does not require cleaning. Keep fluids away from the POC-W22A-H81 interior.
- Never drop any objects or liquids through the openings of the POC-W22A-H81.

### B.2.2 Cleaning Tools

Some components in the POC-W22A-H81 may only be cleaned using a product specifically designed for the purpose. In such case, the product will be explicitly mentioned in the cleaning tips. Below is a list of items to use when cleaning the POC-W22A-H81.

- **Cloth** – Although paper towels or tissues can be used, a soft, clean piece of cloth is recommended when cleaning the POC-W22A-H81.
- **Water/Ethanol alcohol** – A cloth moistened with water or 75% ethanol alcohol can be used to clean the POC-W22A-H81.
- **Using solvents** – The use of solvents is not recommended when cleaning the POC-W22A-H81 as they may damage the plastic parts.
- **Cotton swabs** - Cotton swabs moistened with water are excellent tools for wiping hard to reach areas.
- **Foam swabs** - Whenever possible, it is best to use lint free swabs such as foam swabs for cleaning.

Appendix

C

# EMC Test Summary

---

## POC-W22A-H81 Medical Panel PC

<b>Guidance and manufacturer's declaration – electromagnetic emissions</b>		
<p>The model POC-W22A-H81 is intended for use in the electromagnetic environment specified below. The customer or the user of the model POC-W22A-H81 should assure that it is used in such an environment.</p>		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11		The model POC-W22A-H81 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11		The model POC-W22A-H81 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2		
Voltage fluctuations/ flicker emissions IEC 61000-3-3		

**Recommended separation distances between portable and mobile RF communications equipment and the model POC-W22A-H81**

The model POC-W22A-H81 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model POC-W22A-H81 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model POC-W22A-H81 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## POC-W22A-H81 Medical Panel PC


<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The model POC-W22A-H81 is intended for use in the electromagnetic environment specified below. The customer or the user of the model POC-W22A-H81 should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact  ±8 kV air	±6 kV contact  ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst  IEC 61000-4-4	±2 kV for power supply lines  ±1 kV for input/output lines	±2 kV for power supply lines  ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s)  ±2 kV line(s) to earth	±1 kV line(s) to line(s)  ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
interruptions and voltage variations on power supply input lines  IEC 61000-4-11	<5 % <i>UT</i> (>95 % dip in <i>UT</i> ) for 0,5 cycle  40 % <i>UT</i> (60 % dip in <i>UT</i> ) for 5 cycles  70 % <i>UT</i> (30 % dip in <i>UT</i> ) for 25 cycles  <5 % <i>UT</i> (>95 % dip in <i>UT</i> ) for 5 sec	<5 % <i>UT</i> (>95 % dip in <i>UT</i> ) for 0,5 cycle  40 % <i>UT</i> (60 % dip in <i>UT</i> ) for 5 cycles  70 % <i>UT</i> (30 % dip in <i>UT</i> ) for 25 cycles  <5 % <i>UT</i> (>95 % dip in <i>UT</i> ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model POC-W22A-H81 requires continued operation during power mains interruptions, it is recommended that the model POC-W22A-H81 be powered from an uninterruptible power supply or a battery.



Power frequency (50/60 Hz) magnetic field  IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: <i>UT</i> is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity			
The model POC-W22A-H81 is intended for use in the electromagnetic environment specified below. The customer or the user of the model POC-W22A-H81 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the model POC-W22A-H81, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  <b>Recommended separation distance</b> $d = 1,2\sqrt{P}$  $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz  $d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz  where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in metres (m).
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	V/m	

**POC-W22A-H81 Medical Panel PC**

			<p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
--	--	--	---

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model POC-W22A-H81 is used exceeds the applicable RF compliance level above, the model POC-W22A-H81 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the model POC-W22A-H81.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than V/m.

Appendix

**D**

# BIOS Menu Options

---

## POC-W22A-H81 Medical Panel PC

System Date [xx/xx/xx] .....	79
System Time [xx:xx:xx] .....	79
ACPI Sleep State [S1 only (CPU Stop Clock)] .....	80
Wake System with Fixed Time [Disabled] .....	81
Security Device Support [Enable] .....	83
Hyper-threading [Enabled] .....	84
Active Processor Cores [All] .....	84
Intel Virtualization Technology [Disabled] .....	85
EIST [Enabled] .....	85
SATA Controller(s) [Enabled] .....	86
SATA Mode Selection [IDE] .....	86
USB Devices .....	87
Legacy USB Support [Enabled] .....	87
Serial Port [Enabled] .....	88
Change Settings [Auto] .....	88
Device Mode [RS232] .....	89
Serial Port [Enabled] .....	89
Change Settings [Auto] .....	89
Hardware Health Status .....	90
CPU_FAN1 Smart Fan Control/SYS_FAN1 Smart Fan Control [Auto Mode] .....	91
Fan start/off temperature .....	92
Fan start PWM .....	92
Fan slope PWM .....	92
Console Redirection [Disabled] .....	93
Auto Recovery Function [Disabled] .....	94
Restore AC Power Loss [Last State] .....	96
PCIe Speed [Auto] .....	96
Detect Non-Compliance Device [Disabled] .....	97
Azalia(HD Audio) [Enabled] .....	98
VT-d [Disabled] .....	98
Primary Display [Auto] .....	99
DVMT Pre-Allocated [256M] .....	99
DVMT Total Gfx Mem [MAX] .....	100
Primary IGFX Boot Display [VBIOS Default] .....	100

Bootup NumLock State [On].....	102
Quiet Boot [Enabled] .....	102
Option ROM Messages [Force BIOS].....	102
Launch PXE OpROM [Disabled] .....	103
UEFI Boot [Disabled] .....	103
Administrator Password .....	104
User Password .....	104
Save Changes and Reset .....	104
Discard Changes and Reset .....	104
Restore Defaults .....	105
Save as User Defaults .....	105
Restore User Defaults .....	105
SEL Components [Enabled].....	106
Erase SEL [No] .....	106
When SEL is Full [Do Nothing].....	107
Log EFI Status Codes [Both] .....	107
Configuration Address source [Unspecified] .....	108



Appendix

E

# Watchdog Timer

---



**NOTE:**

The following discussion applies to DOS. Contact IEI support or visit the IEI website for drivers for other operating systems.

The Watchdog Timer is a hardware-based timer that attempts to restart the system when it stops working. The system may stop working because of external EMI or software bugs. The Watchdog Timer ensures that standalone systems like ATMs will automatically attempt to restart in the case of system problems.

A BIOS function call (INT 15H) is used to control the Watchdog Timer.

INT 15H:

<b>AH – 6FH Sub-function:</b>	
AL – 2:	Sets the Watchdog Timer’s period.
BL:	Time-out value (Its unit-second is dependent on the item “Watchdog Timer unit select” in CMOS setup).

**Table E-1: AH-6FH Sub-function**

Call sub-function 2 to set the time-out period of Watchdog Timer first. If the time-out value is not zero, the Watchdog Timer starts counting down. When the timer value reaches zero, the system resets. To ensure that this reset condition does not occur, calling sub-function 2 must periodically refresh the Watchdog Timer. However, the watchdog timer is disabled if the time-out value is set to zero.

A tolerance of at least 10% must be maintained to avoid unknown routines within the operating system (DOS), such as disk I/O that can be very time-consuming.

## POC-W22A-H81 Medical Panel PC

**NOTE:**

The Watchdog Timer is activated through software. The software application that activates the Watchdog Timer must also deactivate it when closed. If the Watchdog Timer is not deactivated, the system will automatically restart after the Timer has finished its countdown.

**EXAMPLE PROGRAM:**

**; INITIAL TIMER PERIOD COUNTER**

;

**W\_LOOP:**

;

```

MOV      AX, 6F02H      ;setting the time-out value
MOV      BL, 30         ;time-out value is 48 seconds
INT      15H

```

;

**; ADD THE APPLICATION PROGRAM HERE**

;

```

CMP      EXIT_AP, 1     ;is the application over?
JNE      W_LOOP        ;No, restart the application

```

```

MOV      AX, 6F02H     ;disable Watchdog Timer
MOV      BL, 0        ;
INT      15H

```

;

**; EXIT ;**

Appendix

**F**

# Hazardous Materials Disclosure

---

## POC-W22A-H81 Medical Panel PC

The details provided in this appendix are to ensure that the product is compliant with the Peoples Republic of China (China) RoHS standards. The table below acknowledges the presences of small quantities of certain materials in the product, and is applicable to China RoHS only.

A label will be placed on each product to indicate the estimated “Environmentally Friendly Use Period” (EFUP). This is an estimate of the number of years that these substances would “not leak out or undergo abrupt change.” This product may contain replaceable sub-assemblies/components which have a shorter EFUP such as batteries and lamps. These components will be separately marked.

Please refer to the following table.

Part Name	Toxic or Hazardous Substances and Elements					
	Lead (Pb)	Mercury (Hg)	Cadmium (Cd)	Hexavalent Chromium (CR(VI))	Polybrominated Biphenyls (PBB)	Polybrominated Diphenyl Ethers (PBDE)
Housing	O	O	O	O	O	O
Display	O	O	O	O	O	O
Printed Circuit Board	O	O	O	O	O	O
Metal Fasteners	O	O	O	O	O	O
Cable Assembly	O	O	O	O	O	O
Fan Assembly	O	O	O	O	O	O
Power Supply Assemblies	O	O	O	O	O	O
Battery	O	O	O	O	O	O
<p>O: This toxic or hazardous substance is contained in all of the homogeneous materials for the part is below the limit requirement in SJ/T11363-2006 (now replaced by GB/T 26572-2011).</p> <p>X: This toxic or hazardous substance is contained in at least one of the homogeneous materials for this part is above the limit requirement in SJ/T11363-2006 (now replaced by GB/T 26572-2011).</p>						



此附件旨在确保本产品符合中国 RoHS 标准。以下表格标示此产品中某有毒物质的含量符合中国 RoHS 标准规定的限量要求。

本产品上会附有“环境友好使用期限”的标签，此期限是估算这些物质“不会有泄漏或突变”的年限。本产品可能包含有较短的环境友好使用期限的可替换元件，像是电池或灯管，这些元件将会单独标示出来。

部件名称	有毒有害物质或元素					
	铅 (Pb)	汞 (Hg)	镉 (Cd)	六价铬 (CR(VI))	多溴联苯 (PBB)	多溴二苯 醚 (PBDE)
壳体	○	○	○	○	○	○
显示	○	○	○	○	○	○
印刷电路板	○	○	○	○	○	○
金属螺帽	○	○	○	○	○	○
电缆组装	○	○	○	○	○	○
风扇组装	○	○	○	○	○	○
电力供应组装	○	○	○	○	○	○
电池	○	○	○	○	○	○

○: 表示该有毒有害物质在该部件所有物质材料中的含量均在 SJ/T 11363-2006 (现由 GB/T 26572-2011 取代) 标准规定的限量要求以下。

X: 表示该有毒有害物质至少在该部件的某一均质材料中的含量超出 SJ/T 11363-2006 (现由 GB/T 26572-2011 取代) 标准规定的限量要求。